


# **EYE INSTRUMENTS COMPANY**



ALBARIARTPRINTER  
0300-4151546

**Straight Tying Forceps**

1401011W  
1401012W



Straight shafts with 4.5mm tying platforms, 85mm long

		<u>Tips</u>
1401011W		Standard
1401012W		Ultrafine

**Mcpherson Tying Forceps**

1401013W  
1401014W  
1401015W  
1401016W



Angled shafts with 4.5mm tying platforms, 85mm long

		<u>Tips</u>	<u>Shafts</u>
1401013W		Standard	12mm
1401014W		Ultrafine	7mm
1401015W		Ultrafine	10mm
1401016W		Ultrafine	12mm

**Castroviejo Tying Forceps**

1401031W  
1401032W




5.0mm long smooth jaws for 8/0 to 11/0 suture. Shafts available straight or angled 45°, Castroviejo handle. Overall length 108mm

		<u>Shafts</u>
1401031W		Straight
1401032W		Angled

**Straight Tying Forceps**

1401021W  
1401022W



Straight shafts with 4.5mm tying platforms, 115mm long

		<u>Tips</u>
1401021W		Standard
1401022W		Ultrafine

**Mcpherson Tying Forceps**

1401023W  
1401024W




Angled 45° shafts with 4.5mm tying platforms 115mm long

		<u>Tips</u>
1401023W		Standard
1401024W		Ultrafine

ALBARIARTPRINTER  
0300-4151546

**Girard Tying Forceps**

1401041W  
1401042W



4.5mm tying platforms.  
Round knurled handle with guide pin,  
105mm long

	<u>Shafts</u>
1401041W	Straight
1401042W	Curved

**Curved Tying Forceps**

1401018W  
1401026W




Curved shafts with 6mm tying platforms.  
Ideal for 10/0 or 11/0 suture.

	<u>Overall Length</u>
1401018W	85mm
1401026W	115mm

**Toothed**

**Straight Toothed Forceps**

1401111W  
1401112W



1 × 2 teeth, with 5.0mm tying platforms  
85mm long

	<u>Teeth</u>
1401111W	0.12mm
1401112W	0.10mm

**Colibri Toothed Forceps**

1401113W  
1401123W




1 × 2 teeth, 0.12mm, with tying platforms,  
85mm long

	<u>Overall Length</u>
1401113W	85mm
1401123W	115mm

**Kelman-Mcpherson Toothed Forceps**

1401114W



1 × 2 teeth, 0.12mm, with 5.0mm tying platforms, and 10mm angled shafts.  
85mm long

ALBARIARTPRINTER  
0300-4151546

**Botvin Iris Toothed Forceps**




1401116W



1 × 2 teeth, 0.10mm for peripheral iridectomies, 85mm long

**Colibri Conjunctiva Toothed Forceps**




1401117W



1 × 2 teeth, 0.12mm, with 5.0mm tying platforms, 85mm long

**Bonn Toothed Forceps**




1401160W



1 × 2 teeth, 0.12mm, 85mm long

**Straight Toothed Forceps**




1401121W  
1401122W



1 × 2 teeth, with 5.0mm tying platforms,  
115mm long

<u>Teeth</u>	
1401121W	0.12mm
1401122W	0.10mm

**Toothed Forceps**



1401131W  
1401132W



1 × 2 teeth, with 5.0mm tying platforms,  
108mm long

<u>Teeth</u>	
1401131W	0.12mm
1401132W	0.40mm

**Harms-Colibri Toothed Forceps**



1401151W  
1401161W



1 × 2 teeth, with tying platforms,  
75mm long

<u>Teeth</u>	
1401151W	0.12mm
1401161W	0.30mm

ALBARIARTPRINTER  
0300-4151546



**Straight Toothed Forceps**

1401152W 1401154W			<b>Overall Length</b>	
			1401152W	85mm
		1 × 2 teeth, 0.12mm with tying platforms	1401154W	115mm

**Colibri Toothed Forceps**

1401153W				
		1 × 2 teeth, 0.12mm, with tying platforms, 85mm long		



**Bishop-Harmon Tissue Forceps**

1401156W				
		1 × 2 teeth, delicate, 85mm long		



**Castroviejo Suturing Forceps**

1401157W 1401158W			<b>Teeth</b>	
			1401157W	0.12mm
		1 × 2 teeth, with 5.0mm tying platforms 108mm long	1401158W	0.30mm



**Straight Toothed Forceps**

1401141W				
		1 × 2 teeth, 0.12mm, with 5.0mm tying platforms, 105mm long		

**Tissue Forceps**

1401159W 1401162W			<b>Shafts</b>	
			1401159W	Straight
		1 × 2 teeth, delicate, 100mm long	1401162W	Curved

**Nucleus Removal Forceps**

1401155W				
		Two rows of teeth on jaws provide a firm purchase on the nucleus fragment for removal, 100mm long		

ALBARIARTPRINTER  
0300-4151546

**Polack Corneal Suturing Forceps**

1401372W



1mm spread, very delicate 1 × 2 teeth, 0.12mm, 115mm long

**Polaok Corneal Suturing Forceps**

1401381W



Colibri style. Two 2.75mm long tips, separated 1mm, with 1 × 2 teeth, 75mm long

**Notched**

**Pierse-Colibri Corneal Forceps**


1401101W



0.25mm notched, with tying platforms, 75mm long

**Straight Notched Forceps**

1401061W  
1401062W




0.25mm notched, with tying platforms,  
85mm long

**Tips**

1401061W	Standard
1401062W	Ultrafine

**Angled Notched Forceps**

1401063W



0.25mm notched, with tying platforms, 85mm long

**Curved Notched Forceps**

1401065W




0.25mm notched, with tying platforms, 85mm long

ALBARIARTPRINTER  
0300-4151546

**Straight Notched Forceps**

1401071W  
1401072W




0.25mm notched, 115mm long

Tips	
1401071W	Standard
1401072W	Ultrafine

**Pierse Corneal Forceps**

1401081W  
1401082W



0.25mm notched, very delicate  
115mm long

Shafts	
1401081W	Straight
1401082W	Curved

**Hofmann-Polack Corneal Suture Forceps**

1401371W



1mm spread, 0.25mm notched, 115mm long

**Capsulorhexis**

**Bloomberg Pierse Capsulorhexis Forceps**

1401171W  
1401189W




45° angled straight shafts.  
Tip to angled length 14.0mm  
used to perform continuous tear  
circular or oval capsulotomies.

Overall Length	
1401171W	85mm
1401189W	115mm

**Bloomberg Pierse Capsulorhexis Forceps**

1401172W  
1401188W




45° angled gently curved shafts.  
Tip to angled length 14.0mm  
used to perform continuous tear  
circular or oval capsulotomies.

Overall Length	
1401172W	85mm
1401188W	115mm

ALBARTPRINTING  
0300-151546

**Ultrata Style Capsulorhexis Forceps**

1401173W  
1401174W




Fine angled sharp tips, 12mm 45° angled shafts, 85mm long

	<u>Shafts</u>
<b>1401173W</b>	Straight
<b>1401174W</b>	Curved

**Bloornberg Pierson Capsulorhexis Forceps**

1401175W  
1401187W




Capsulorhexis notched.  
45° angled tips with curved 12mm shafts, with an iris stop

	<u>Overall Length</u>
<b>1401175W</b>	85mm
<b>1401187W</b>	115mm

**Corydon Style Capsulorhexis Forceps**

1401176W  
1401186W



Fine 45° angled sharp tips with curved 12mm shafts for performing capsulorhexis through tunnel incision. Iris stops to prevent inadvertent iris capture. Cystotome teeth can be used to create initial flap in anterior capsule for single instrument capsulorhexis procedure

	<u>Overall Length</u>
<b>1401176W</b>	85mm
<b>1401186W</b>	115mm

**Lehner-Ultrata Style Capsulorhexis Forceps**

1401221W



Fine angled sharp tips with curved 12mm shafts for performing capsulorhexis through tunnel incisions. With guide pin, 105mm long.

ALBANYARTPRINTER  
0300-4151546



**Bloomberg Pierson Capsulorhexis Forceps**

1401222W



Used to perform continuous tear circular or oval capsulotomies.  
 12mm 45° angled curved shafts with an iris stop.  
 With guide pin, 105mm long

**Corydon Capsulorhexis Forceps**

1401223W



Round handle for finger tip rotational control. Fine 45° angled tips with curved 12mm shafts and iris stop. Cystotome teeth can be used to create initial flap in anterior capsule for single instrument capsulorhexis procedure. With guide pin, 105mm long

**O'Gawa-Utrata Capsulorhexis Forceps**

1401224W



Fine 45° angled tips with shafts 12mm with an iris stop. Round handle for finger tip rotational control. Cystotome teeth be used to create initial flap in anterior capsule for single instrument capsulorhexis procedure. With guide pin, 105mm long

**Capsulorhexis Forceps**

1401225W



Cystotome tips can be use to create initial flap in anterior capsule for single instrument capsulorhexis procedure. 45° angled 12mm straight shafts and Iris stop, 120mm long

**Capsulorhexis Forceps**

1401226W



Cystotome tips can be used to create initial flap in anterior capsule for single instrument capsulorhexis procedure. 45° angled 12mm curved shafts and Iris stop for performing capsulorhexis through tunnel incisions, 120mm long

ALBARART-PRINTER  
 0300-4151546

**Utrata Style Capsulorhexis Forceps**

1401227W



Round handle for finger tip rotational control. Iris stop to prevent inadvertant iris capture. 45° angled tips with slim 10mm shafts to fit through phaco incision, without wound tenting. With guide pin, 105mm long

**Nevyas Capsulorhexis Forceps**

1401201W



Fine 45° angled tips with 11mm straight shafts and iris stop. Cystotome teeth can be used to create initial flap in anterior capsule for single instrument capsulorxis procedure. 108mm long

**IOL**

**Kelman-Mcpherson Forceps**

1401261W

1401271W



Angled 45° shafts. With 7.5mm long smooth jaws.

1401261W

1401271W

**Overall Length**

85mm

115mm

**Kelman-Mcpherson Lens Forceps**

1401301W



7.5mm long smooth jaws. Very delicate, 108mm long

**Bechert Lens Forceps**

1401302W




10mm long smooth jaws. Very delicate, 108mm long

ALBARIARTPRINTER  
0300-4151546

**Clayman Lens Forceps**

1401281



Gently curved shafts with 90° angled 1.0mm tips. 120mm long

**Lens Forceps**

1401282



Curved shafts with 1mm ring tips, 120mm long

**Kratz Style Lens Forceps**

1401283



Curved narrow smooth jaws. Very delicate, 120mm long.

**Jaffe-Blaydes Lens Forceps**

1401284



Angled 30° narrow smooth jaws. Very delicate, 115mm long.

**Blaydes Angled Lens Forceps**

1401285



8mm 45° angled shafts with smooth jaws, 115mm long

**Shepard Lens Forceps**

1401286



Gently curved shafts with flat tungsten carbide coated jaws, 125mm long

ALBARIARTPRINTER  
0300-4151546

**Anis Lens Forceps**

1401287



9mm long smooth jaws, angled 45° ,  
extra delicate, 115mm long

**Kelman-Mcpherson Lens Forceps**

1401321



Lens holding Forceps. Smooth jaws for holding soft IOLs,  
102mm long

**Lens Implantation Forceps**

1401331



The jaw's smooth surfaces and rounded edges grip  
the lens firmly yet gently, 115mm long

**Lens Implantation Forceps**

1401332




Reverse cross-action. Smooth rounded curved jaws, 105mm long

**Lens Implantation Forceps**

1401333

1401334



Standard cross-action. Smooth and curved jaws,  
108mm long

**1401333**

**1401334**


**Lock**

Without lock

With lock

**Lens Implantation Forceps**

1401335



Standard cross-action. The thin low-profile jaws securely grip the optic  
and avoid lens migration. With lock, 108mm long

ALBARIARTPRINTER  
03004151546

**Lens Implantation Forceps**


1401336



The ultra-fine, smooth and well rounded blades can grasp the optic low down without lens migration. Thin profile blades with convex inner surfaces, 108mm long

**Seibel Style Lens Folding Forceps**

1401311  
1401312



Folding 6mm IOL's. Paddle style jaws with handle stop pin to prevent over folding IOL, 115mm long

**1401311**  
**1401312**

**Lock**  
With lock  
Without lock

**Seibel Style Lens Folding Forceps**

1401313



Incorporates tying platform tips for tucking the haptics into folded 6 mm IOL. Handle stop pin to prevent over folding IOL, 115mm long

**Seibel Style Lens Folding Forceps**


1401314



Folding 5.5/6.0mm IOL's. Handle stop pin to prevent over folding IOL, 115mm long

**Steinhart Style Lens Folding Forceps**

1401315



Triangular shaped inner jaws, handle stop pin to prevent over folding IOL, 115mm long

**Lens Folding Forceps**

1401327




Universal lens folding forceps, for 5.5/6.0mm IOL's. The box design of the jaws folds and encapsulates the IOL and the side cut outs enable the perfect transfer to the inserter, 108mm long

ALBARIART PRINTER  
0300-4151546

**Lens Folding Forceps**

1401328



Design facilitates pick up and folding of 5.5/6.0mm IOL's from a flat surface and incorporates a handle stop pin to prevent over folding, 108mm long


**Lens Folding Forceps**

1401329



Universal lens folding forceps, for 5.5/6.0mm IOL's. The box design of the jaws folds and encapsulates the IOL and the side cut outs enable the perfect transfer to the inserter, 108mm long

1401338W **Lens Implantation Forceps**



**Prechopper**

**Akahoshi Prechopper**



1401391



Universal smooth pointed tips, 12mm angled shafts, cross-action, maximum tip opening 3.0mm, 120mm long

**Akahoshi Prechopper**

1401392



Universal smooth pointed tips, with 3mm cutter, 12mm angled shafts, cross-action, maximum tip opening 2.5mm, 120mm long

**Akahoshi Prechopper**

1401397



Conventional smooth sharp pointed tips, 12mm angled shafts cross-action, maximum tip opening 3.0mm, 120mm long

**Akahoshi Prechopper**

1401398



Conventional smooth sharp pointed tips. Straight shafts cross-action, maximum tip opening 3.0mm, 120mm long

ALBARIARTPRINTER  
0300-4151546

**Akahoshi Prechopper**



1401399



Conventional smooth sharp pointed tips, 12mm angled curved shafts cross-action, maximum tip opening 3.0mm, 120mm long

**Akahoshi Phaco Prechopper**



1401402

1401403



4.5mm long delicate jaws with conventional smooth sharp tips for cracking dividing and rotating the nucleus without sculpting, 120mm long

<b>Shafts</b>	
<b>1401402</b>	Straight
<b>1401403</b>	Angled

**Nucleus Bisector/Vectis Forceps**



1401393

**Nucleus Trisector/Vectis Forceps**




1401394

110mm long

**Conjunctiva**

**Fechtner Conjunctiva Forceps**



1401400



Very delicate ring tips, with tying platforms, 112mm long

**Chalazion**

**Chalazion Forceps**



1401405

Oval 12 × 14mm, 98mm long



1401406

1401407

1401408



**1401407** Round 15mm, 98mm long



**1401408** Oval 17 × 9mm, 93mm long

ALBARIART-PRINTER  
0300-4151546

**Chalazion Forceps**

1401409



1401412 **Cilia Forceps**





72mm long

1401416 **Scleral Plug Forceps**



Designed to securely hold scleral plugs during insertion and removal.  
The cross-action handle reduces hand fatigue while grasping plug.  
118mm long

1401414 **Hartmann Straight Hemostatic Mosquito Forceps**



Straight 20mm serrated jaws, ring handle with ratchet lock,  
95mm long

1401415 **Hartmann Curved Hemostatic Mosquito Forceps**



Curved 20mm serrated jaws, 95mm long

ALBANYARTPRINTER  
0300-4151546



1401417 Baby Jones Towel Clamps



65mm long

Landers Vitrectomy Lens Forceps

1401418



For use with the Landers Vitrectomy lens set. Shafts angled 45° with curved grooved blades for holding the Landers Vitrectomy Lens, 115mm long

Jewelers

1401421 Jeweler Style Forceps



0.1 × 0.06mm tips, 110mm long

1401422 Jeweler Style Forceps



0.05 × 0.01mm tips, 90° bent, 110mm long

1401423 Jeweler Style Forceps



0.05 × 0.01mm tips, bent at 45°, 110mm long

1401424 Jeweler Style Forceps





0.07 × 0.03mm tips, 110mm long

ALBARIARTPRINTER  
0300-4151546

1401425 Lens Folding Forceps




1401611 Toothed Forceps




1401613 Curved Tying Forceps



1401614 Straight Tying Forceps



1401616 Tissue Forceps



1401641 Shepard Lens Forceps



1401642 Shepard Lens Forceps



ALBARIARTPRINTER  
0300-4151546

# Needle Holders

## Curved Needle Holder

1402011W  
1402016W



Curved, 11mm medium jaws, sharp tips,  
115mm long

1402011W  
1402016W

Lock
Without lock
With lock



## Barraquer Needle Holder-Curved

1402101W  
1402106W



Delicate, 9mm gently curved jaws, sharp tips,  
115mm long

1402101W  
1402106W

Lock
Without lock
With lock



## Troutman Needle Holder-Curved

1402231W  
1402232W



10mm gently curved jaws, blunt tips,  
115mm long

1402231W  
1402232W

Lock
Without lock
With lock



ALBARIARTPRINTER  
0300-4151546

# Needle Holders

## Curved Needle Holder

1402233W  
1402241W




Curved, 6mm delicate jaws, 110mm long

1402233W  
1402241W

**Lock**  
Without lock  
With lock

## Straight Needle Holder

1402234W  
1402242W



Straight, 6mm delicate jaws, 110mm long

1402234W  
1402242W

**Lock**  
Without lock  
With lock

## Bechart-Sinskey Needle Holder

1402235W  
1402243W



Delicate, 9mm gently curved jaws, sharp tips,  
115mm long

1402235W  
1402243W

**Lock**  
Without lock  
With lock

## Side Curved Needle Holder

1402236W



Blunt tips,  
115mm long

## Fine Curved Needle Holder

1402237W  
1402244W



Curved, 9mm fine jaws, sharp tips,  
115mm long

1402237W  
1402244W

**Lock**  
Without lock  
With lock

ALBARIARTPRINTER  
0300 451546

# Needle Holders

## Barraquer Needle Holder

1402238W  
1402245W



Very delicate curved jaws, sharp tips,  
115mm long

1402238W  
1402245W

Lock
With lock
Without lock



## Azar Needle Holder

1402176W  
1402177W



Fine, 9mm gently curved jaws,  
115mm long

1402176W  
1402177W

Lock
Without lock
With lock




## Castroviejo Needle Holder

1402191W  
1402192W

Gently curved delicate jaws, blunt tips,  
140mm long

1402191W  
1402192W

Lock
Without lock
With lock



## Castroviejo Needle Holder

1402193W  
1402195W

Standard, 11mm curved jaws, blunt tips,  
140mm long

1402193W  
1402195W

Lock
With lock
Without lock




ALBARIARTPRINTER  
0300-4151546

# Needle Holders

## Castroviejo Needle Holder

1402204W  
1402205W



Delicate jaws, Straight blunt tips,  
140mm long

1402204W  
1402205W

Lock	
1402204W	Without lock
1402205W	With lock

## Castroviejo Needle Holder

1402203W  
1402206W



Standard, 11mm jaws, blunt tips,  
140mm long

1402203W  
1402206W

Lock	
1402203W	Without lock
1402206W	With lock

## Castroviejo Needle Holder

1402207W



## Castroviejo Needle Holder


1402211W  
1402212W



Fine 9mm curved jaws,  
wide serrated handle, sharp tips,  
125mm long

1402211W  
1402212W


Lock	
1402211W	Without lock
1402212W	With lock



ALBARRARTPRINTER  
0300-4151546

## Castroviejo Needle Holder

1402213W  
1402214W



Fine 9mm straight jaws,  
wide serrated handle, sharp tips,  
125mm long

### Lock

1402213W Without lock  
1402214W With lock



## Castroviejo Needle Holder

1402405



130mm long

## Castroviejo Needle Holder

1402415



115mm long


ALBARIARTPRINTER  
0300-4151546

## 1403026 Mcpherson-Westcott Conjunctiva Scissors



Curved, blunt tips.  
Tip to pivot length 10.5mm,  
105mm long

## 1403071R Conjunctiva Scissors



Gently curved blades. Tip to pivot length 12mm. Blunt tips,  
115mm long

## 1403056 Osher Universal Scissors



Thin curved blades. Blunt tips, 21mm from pivot to tip, 115mm long

## Corneal

## Castroviejo Corneal Section Scissors

1403027



Right


1403028



Left

Blunt tips and lower blade  
0.5mm longer than the upper blade, 13mm from pivot to tip, 105mm long

## 1403029 Castroviejo Keratoplasty Scissors



Universal with blunt angled blades. Blunt tips and lower blade longer than  
the upper blade, 16mm from pivot to tip, 105mm long



**Castroviejo Corneal Section Scissors**

1403043




Right

1403044




Left

Blunt tips and lower blade 0.5mm longer than the upper blade,  
9mm from pivot to tip, 105mm long

**1403030R Castroviejo Curved Corneal Scissors**


Universal medium model with blunt tips. 16mm from pivot to tip,  
116mm long

**1403031 Castroviejo Corneal Scissors**

Curved blunt tips, small blades, 105mm long


**1403072R Castroviejo Corneal Scissors**

16mm from pivot to tip, 116mm long

**1403073R Castroviejo Corneal Scissors**


Blunt tips, curved blades, 11.0mm from pivot to tip, 112mm long

1403075 **Castroviejo Curved Scissors**



Blunt tips, tip to pivot length 11.0mm, 105mm long

1403076 **Castroviejo Straight Scissors**



Blunt tips, tip to pivot length 11.0mm, 105mm long

Iris

1403032 **Iris Scissors**



Sharp tips, straight blades, cutting length 4.0mm, blade length 11.0mm, 90mm long

1403033 **Iris Scissors**



Sharp tips, curved blades, cutting length 4.0mm, blade length 11.0mm, 90mm long

1403034 **Iris Scissors**



Sharp tips, curved blades, blade length 11.0mm, 105mm long

**Iris Scissors**

1403045



Sharp tips, 11.0mm blades, 105mm long

ALBARIARTPRINTER  
0300-4151546

1403008 Mcpherson-Vannas Curved Iris Scissors



Stainless steel. Sharp pointed tips, 12mm from pivot to tip, 90mm long

1403077 Barraquer Iris Scissors



Stainless steel. Blunt tips, 8mm from bend to tip. Squeeze action handle. Overall length 60mm.


Stitch

1403059 Mcpherson-Westcott Stitch Scissors



Very sharp pointed tips, 10.5mm from pivot to tip, 105mm long

1403078 Westcott Stitch Scissors



Very sharp pointed tips, 19mm from pivot to tip, overall length 115mm

Mcpherson-Westcott Stitch Scissors

1403074R  
1403091




Sharp pointed tips,  
13mm from pivot to tip

1403074R  
1403091

Overall Length	
1403074R	115mm
1403091	140mm

Tenotomy


1403058 Westcott Tenotomy Scissors



Curved blunt tips, 21mm from pivot to tip, 115mm long

ALBARIARTPRINTER  
0300-4151546

1403079 **Westcott Tenotomy Scissors**




Curved blunt tips, 19mm from pivot to tip, 110mm long

Utility

**Eye Scissors**

1403062



1403064



Pointed tips, 32mm from pivot to tip, 110mm long

	<u>Blades</u>
1403062	Straight
1403064	Curved

Capsulotomy

**Vannas Capsulotomy Scissors**

1403035

1403036




Sharp pointed tips, straight blades, 108mm long

	<u>Blades</u>
1403035	10.5mm blades
1403036	12.5mm blades

**Vannas Capsulotomy Scissors**

1403037

1403038



Sharp pointed tips, curved blades, 95mm long

	<u>Blades</u>
1403037	10.5mm blades
1403038	12.5mm blades

ALBARIARTPRINTER  
03002451546

Vannas Capsulotomy Scissors

1403039  
1403040




Sharp pointed tips. Angled blades, 95mm long

**Blades**

1403039	10.5mm blades
1403040	12.5mm blades

Vannas

1403001 Vannas Scissors



Sharp pointed tips. Straight 7.5mm blades, 9.5mm from pivot to tip, 82mm long

1403002 Vannas Scissors



Sharp pointed tips. Curved 7.5mm blades, 9.5mm from pivot to tip, 82mm long

1403007 Vannas Scissors



Sharp tips. Angled to side 7.0mm blades, 9.0mm from pivot to tip, 82mm long

1403061 Vannas Scissors



Extra thin 7mm long blades. Sharp points, 9.0mm from pivot to tip, 105mm long

Vannas Scissors

1403041



1403042



Extra thin blades, sharp points, from pivot to tip 15mm, 108mm long


**Blades**

1403041	Straight
1403042	Curved

ALBARIARTPRINTER  
030-4151546

**Kuglen Iris Hook**

1404101  
1404102



Iris hook and lens manipulator



Shaft	
1404101	Angled
1404102	Straight

**Jaffe-Knolle Iris Hook**

1404108  
1404109



Iris hook, very blunt, small



Shaft	
1404108	Angled
1404109	Straight

**Fukasaku Small Pupil Snapper Hook**

1404143



With a smoothly polished notch superiorly to push back the iris with sharp inner edge for splitting the nucleus



**IOL Hooks**

**Fenzel Lens Manipulating Hook**

1404103  
1404104



0.15mm dia. T-shaped tip



Shaft	
1404103	Angled
1404104	Straight

**Jaffe-Bechert Lens Nucleus Rotator**

1404105



Blunt forked tip



**Lester IOL Manipulator**

1404106  
1404107



IOL manipulator for lens with or without holes



Shaft	
1404106	Angled and curved
1404107	Straight and curved

**Sinskey Lens Manipulating Hook**

1404110





0.15mm dia blunt tip, angled to right





ALBARIARTPRINTER  
0300-4151546


**Lewicky Lens Manipulating Hook**

1404111    
 0.15mm diameter blunt tip, angled to right on a curved shaft



**Rentsch Lens Manipulator**

1404112    
 Iris/capsule retractor or pusher/manipulator for IOLs with and without holes

**Maltzman-Fenzl Lens Manipulator**

1404144    
 1404145 0.15mm dia. V-shaped tip **Shaft**  
1404144 Straight  
1404145 Angled

**Repositor**



1404114    
 1404115 Synechia spatula. Semi-sharp tips **Tips**  
1404114 0.7mm  
1404115 0.5mm

**Nucleus Manipulator**



**Drysdale Nucleus Manipulator**

1404113    
 Paddle shaped tip

**Botton Manipulator**

1404135    
 Designed for iris retraction and capsule polishing. With roughened posterior surface, angled 45° shaft, 10mm from tip to bend

**Intracapsular Manipulator**

1404120    
 Designed to safely spin, chop and backcrack a nucleus using less phaco time. Perfect for dialing intraocular lenses in and out of the bag, the wand is also used to provide counter traction in the paracentesis for improved fixation during construction of the phaco incision, capsulorhexis and lens delivery. Designed for manipulating and sweeping nucleus fragments toward a phaco tip in the central safety zone for safe emulsification with less stress to the zonules

ALPHART-PRINTER  
0820-441546

**Akahoshi Nucleus Stainer/manipulator**

1404140



Designed to engage the nucleus equator to stabilize it and provide counter pressure while the prechopper forceps is being inserted into the nucleus core

**Intracapsular Manipulator**

1404150



Designed to safely spin, chop and backcrack a nucleus using less phaco time. Perfect for dialing intraocular lenses in and out of the bag, the wand is also used to provide counter traction in the paracentesis for improved fixation during construction of the phaco incision, capsulorhexis and lens delivery. Designed for manipulating and sweeping nucleus fragments toward a phaco tip in the central safety zone for safe emulsification, with less stress to the zonules

**Intracapsular Manipulator**

1404151



Designed to safely spin, chop and backcrack a nucleus using less phaco time. Perfect for dialing intraocular lenses in and out of the bag, the wand is also used to provide counter traction in the paracentesis for improved fixation during construction of the phaco incision, capsulorhexis and lens delivery. Designed for manipulating and sweeping nucleus fragments toward a phaco tip in the central safety zone for safe emulsification, with less stress to the zonules

**Retina Hooks and Spatulae****Retina Hook**

1404133a



Triangular 0.5mm tip, angled 130°

**Membrane Spatula**


1404133b



Rounded tip spatula

**Retina Spatula**

1404133c



Tip angled 110°

**Retina Hook**

1404133d





# Muscle Hooks

## 1404119 Retinal Detachment Hook



With 1.5mm oval hole. Flattened 13mm

## Helveston Muscle Hook


- 1404152 8mm
- 1404153 10mm
- 1404154 12mm



Delicate muscle hook with blunt spatulated tip

## Helveston Finder Hook

- 1404155 Small
- 1404156 Large



Lightly angled cone shaped tip

## 1404157 Helveston Teaser Hook



6mm hook

## Helveston Retractors

- 1404158 7mm
- 1404159 9mm
- 1404160 11mm



## Helveston Scleral Marking Ruler

1404161



ALBARIARTPRINTER  
0300-4151546

**1404118 Kansas Lens Loop**



3mm wide. Serrated. For small incision nucleus fragment removal

**Kansas Nucleus Bisector**

1404116



For manual small incision surgery

**Kansas Nucleus Trisector**

1404117



For manual small incision surgery

**1404301 Kansas Nucleus Vectis**



Solid used as a "cutting block" for Kansas nucleus bisector and trisector, for left and right handed surgeon

**Phaco Chopper**

1404121

1404122



**Shaft**

1404121 Straight

1404122 Angled

**Nealiann Nucleus Divider**

1404131



Ambidextrous nucleus chopper

**Nichamin Phaco Chopper**

1404123

1404124



Ambidextrous nucleus chopper, tip 1.1mm long.  
Rounded and highly polished distal tip

**Blade**

1404123 Straight

1404124 Curved

ALBARIARTPRINTER  
03004151546

# Choppers

## Phaco Chopper

1404125  
1404125L



Phaco side cutter chopper, rounded and highly polished distal tip

1404125 Right hand  
1404125L Left hand

## Phaco Chopper

1404126  
1404126L



Phaco side cutter chopper, rounded and highly polished distal tip

1404126 Right hand  
1404126L Left hand

## Minami "m" Hook

1404127



Chopper, Splitter, Dissector, Spatula, Rotator and Manipulator in one

## Gillum Nucleus Splitter

1404128  
1404129



1404128

1404129

## Claw Phaco Chopper

1404130



Special triple inner blades construction allows chopping in three directions. Special formed shaft, special rounded and highly polished distal tip

## Phaco Chopper

1404136  
1404137



Ambidextrous nucleus chopper, tip 1.8mm long  
Rounded and highly polished distal tip

	Blade
1404136	Straight
1404137	Curved

## Phaco Chopper

1404138  
1404138L



Phaco side cutter chopper, rounded and highly polished distal tip  
Tip 1.8mm long

1404138 Right hand  
1404138L Left hand

# Choppers

## Haefliger Phaco Cleaver

1404139



Chopper and manipulator in one

## Phaco Chopper

1404141



Special inner blade is offset 45°, rounded distal tip is highly polished, unique elongated distal segment

## Phaco Chopper

1404142

1404142L



Rounded and highly polished distal tip.  
Tip 1.8mm long

1404142 Right hand

1404142L Left hand

## Phaco Chopper

1404302




Paddle manipulator / 60° Phaco chopper



## Phaco Chopper

1404303



With Lieberman Micro finger for nucleus manipulation and phaco chopping.  
With a defined inferior edge for efficient dividing of soft nucleus.  
The "Quick chop" has a sharp point to easily penetrate and split a hard nucleus.  
Double ended for vertical and horizontal phaco chopping



## Schocket Double Ended Scleral Depressor

1404304



Double-end depressor and marker. Round end designed for depressing sclera and curved marking end for reaching behind the globe.  
6mm wide bar and 2 5mm tear drop shape tips



## Phaco Chopper

1404305




1.8mm long



ALBARIARTPRINTER  
0300-4151546

## Helveston Fixation Forceps

1401382



For fixating and retracting the globe during strabismus procedures.

1 x 2 0.5mm teeth. Thumb activated catch for locking.

The instrument features a unique, weighted, curved handle for stabilizing the forceps after affixing it to the globe. A suture may be placed through a hole in the end of the handle and attached to the drape for additional traction, 118mm long

## Fixation Forceps

1401373

1401375



1 x 2 0.12mm teeth at each tip, 115mm long



1401373 13.0mm spread



1401375 3.0mm spread

1404146

## Hofman-Thornton Fixation Ring



Swivel head, 16mm dia

## Fine/Thornton Fixation Ring

1404147

1404148



Swivel head

1404147 13mm/11mm cut out

1404148 16mm/14mm cut out

## Fine Fixation Ring

1404149



Swivel head, 13mm/9mm cut out.

With blunt teeth on both sides

**1405029 Liberman Solid Speculum**



Allow for the protection of the eye lashes and keeps the secretion of the meibomian glands from the operating field, 65mm long

**1405030 Machat Adjustable Speculum**



Provides maximum exposure of the globe for superior LASIK procedure. Open blade design avoids interference of the speculum with the microkeratome, 65mm long

**1405033 Steinert LASIK Speculum LASIK**



Strong with ideal curvature for accommodating the microkeratome suction ring, 65mm long

**1405034 Genisi Speculum LASIK**



Adjustable, strong with ideal curvature for accommodating microkeratome suction ring, 65mm long

**1405027 Cornesi Epithelium Spatula**



**LASIK Flap Manipulator**

1405045



Features a delicate Olive-shaped tip which is used to gently depress the corneal to locate the flap edge, once located engage the edge, with the smoothly polished tip and slides it around the rim to partially open the flap

ALPHA PART PRINTER  
0360-4151546

## Gimbel Stabilization Ring

1405032



Offset 13mm diameter ring with 9mm cut out to accommodate the flap hinge.  
Blunt teeth on both sides for stabilizing either the right or left eye during LASIK

## Moretsky LASIK Hinge Protector Fixation Ring


1405040



The protector is designed to fixate the globe during laser ablation while covering the hinge and flap for complete protection, 11mm inner diameter offers secure fixation of the cornea during ablation without interference. Unique angled protector shields the flap and hinge

## Flap Elevator

1405022



1405025 **Flat Cover**



## Buratto Flap Protector


1405031



Ideal for retracting the flap for maximum exposure of the bed  
and for protecting nasally hinged corneal flap during the laser ablation

## Katzen LASIK Flap Unzipper

1405046



Combines a modified Sinskey hook and a bi-convex spatula with beveled notches for LASIK enhancements

## Caro/Avalos Flap Smoother

1405026



## Flap Smoother

1405021



Designed to conform the corneal flap to the new shape of the stroma bed after laser ablation

## LASIK Multipurpose Forceps

1405028




JAWS  
CROSS  
SECTION

Curved jaws to conform to the curvature of the corneal with spatulated tips and edges to enter, dissect elevate and smooth the flap. It offers firm spring action for dissecting during enhancements and smooth inside jaws surfaces for grasping the flap, This forceps combines the features of a spatula, dissector and forceps, 115mm long

## Perone LASIK Flap Forceps

1405035



Enables delicate lifting of the flap without causing any trauma, 2mm rings, 108mm long

## Manche LASIK Forceps

1405036



Solid oval tips lift flaps for primary and retreatment procedures with minimal risk of tearing or trauma, 108mm long

## Hersh LASIK Retreatment Forceps

1405037



Forceps are designed to lift the LASIK flap easily for retreatment avoiding damage to the flap and facilitating a controlled "peel" of the flap from the stroma. Double grip design to avoid bending the flap and facilitate lifting, 108mm long

ALBANYARTPRINTER  
0300-4151546



## Spatulated LASIK Cannula

1405038



This 26G curved cannula has an end opening for irrigating the interface. Flattened spatulated tip to easily locate and elevate the flap edge.

## Gimbel LASIK Fountain Cannula

1405039



Smooth bullet shaped tip single centerport for directing fluid toward the superior hinge, 25G

## Buratto LASIK Cannula

1405041

1405042



One front and two side ports. Designed for equal flow in three directions

<b>1405041</b>	25G
<b>1405042</b>	23G

1405043

## Banaji LASIK Cannula



Bullet shaped tip. 6 ports for multi-directional irrigation, 25G

1405044


## Vidaurri LASIK Cannula



Double-armed 8 irrigating ports  
25 G

# Speculums

**1405001 Wire Speculum**



With 13.5mm long blades,  
52mm long

**1405002 Lieberman Speculum, V-wire**




The V-shaped wire blades are  
designed to help the surgeon  
retain the drape under the eye lids.  
65mm long

**1405003 Slade LASIK Wire Speculum**



Specifically designed for maximum  
exposure of the globe and to accomodate  
suction ring for efficient operation  
of the keratome, 65mm long

**1405004 The Lindstrom-chu Aspirating Speculum**



Designed to prevent pooling of excess balanced salt solutin from eye  
lid margins through continuous aspiration. Unique reversible speculum  
allowing the surgeon to pivot the device to rest nasally or temporally, 65mm long

**1405005 Lieberman Speculum, K-wire**




Adjustable wire speculum  
with 14.5mm open blades.  
With shaft angled 45°, 65mm long

ALBARTARTPRINTER  
0300-4151546


# Speculums

**1405006 Reversible Lid Speculum**




This speculum can be placed temporarily or nasally, accomodating the suction ring of a mircrokeratome during LASIK surgery. The solid blades keep eye lashes away from the surgical field. While their non-reflective finish eliminates glare from the microscope light, 65mm long

**1405007 Kratz Barraquer Wire Speculum**



With 13.5mm open blades, 52mm long

**1405008 Barraquer Solid Speculum**




With 14.5mm long blades, 52mm long

**1405009 Naqahara Wire Speculum**



Specially designed for temporal approaching cataract surgery. With 13.5mm long blades, 52mm long

**1405010 Grandon Wire Speculum**



It features gently curved back tips to atraumatically retain redundant lid tissue. The heavy gauge coiled fulcrum provides strong spring tension to completely retract and control the patients lids and control redundant lid tissue, 52mm long

ALBARA REPRINTER  
0300-151516


# Speculums

1405011 Glnsberg Eye speculum




Non-reflective blades. With flared ends to control lid skin folds. The lower blade is shorter than the up blade to facilitate easy insertion and removal. The blades are fenestrated to minimize pressure exerted on the globe strongly angled shafts allow the opening and closing mechanism to lie flat along the patient's temple and completely out of the way, 65mm long

1405012 Shapira Speculum




Adjustable, wire opening for phaco and micro lamellar keratoplasty, 65mm long

1405013 Slade LASIK Speculum



52mm long

1405014 Eye Speculum




Swivel open blades to accomodate the variability of patient's anatomical features, 65mm long

1405015 Wire Speculum



65mm long

Squeeze handle  
Instruments for  
Vitreo-Retinal  
surgery



20G, a selection of  
forceps tip are available for Vitreo-Retinal  
procedures as well as anterior segment surgeries,  
such as sphincterotomy, Iridectomy and capsulorhexis.



**1401501**  
Side gripping  
Forceps smooth  
round jaws 2.4mm



**1401502**  
Foreign body  
Forceps with cap  
jaws, straight



Capsulorhexis  
Forceps, 2.5mm  
jaws, 0.7mm  
platforms  
**1401503** Curved  
**1401504** Straight



**1401505**  
Gripping Forceps,  
2.5mm straight



**1401506**  
Gripping Forceps,  
2.0mm smooth  
rounded tips



**1401507**  
Serrated  
straight



**1401508**  
Straight Forceps,  
with 1 x 2 teeth



**1401509**  
Straight, and  
gripping jaws



**1401510W**  
Membrane  
peeler,  
smooth, 45°




**1401511**  
Membrane  
peeler,  
grooved, 45°



**1401512**  
Foreign Body  
Forceps,  
three claws

ALBARIARTPRINTER  
0300-4151546

These squeeze handle instruments feature titanium handles to reduce overall weight. 20G shaft and front ends are made of hardened stainless steel for superior sharpness



**1403151**  
Straight pointed tips  
直尖头剪



**1403152**  
Straight blunt tips  
直钝头剪



**1403153**  
Curved(R12mm)



**1403154**  
Angeld, 45°



**1403155**  
Vertical 40°



**1403156**  
Vertical  
80°



**1403157**  
Vertical, lower  
blade with bead tip

1410001 Blade Breaker



1410002 Blade Handle



1410003 Braunstein Fixed Caliper

1410004



Fine tips mark sclerotomy sites for vitreous surgery

Aphakic/Phakic

1410003	3.5/4.0mm
1410004	3.0/3.5mm

1410005 Castroviejo Caliper



Caliper measures from 0mm to 20mm in 0.5mm increments with scale reading on one side.

## Luntz-Dodick Trabeculectomy Punch

1410011  
1410015  
1410016




With bullet shaped tip for easy insertion into the sclero-corneal tunnel, it's angled cutting edge is designed to engage the scleral lip at the base of the tunnel, 0.5mm deep bite

Diameter Head

1410011	1.0mm
1410015	1.5mm
1410016	1.2mm

## Descemet Membrane Punch Descemet

1410012




1.0mm diameter head, 0.75mm deep bite

ALBARIARTPRINTER  
0300-4151546

## Holth corneoscleral Punch


1410013



1.0mm wide by 7mm long punch

## Christen Punctum Punch

1410014



The punch features a thin tapered tip, which can be inserted into the punctum to excise a 1mm wide section of tissue in a single action. There is no need for multiple cuts, as with a scissors, there by minimizing the amount of bleeding

## Trabeculotome set

1404134



## Urrets-Zavalía sclera marker

1404162



Pressure and mark surface of sclera, blunt tips

## 1410102 Foldable IOL Implantation Holdetr




φ 1.8mm

ALBARIARTPRINTER  
0300-4751546




- 1410301 **Simcoe Double-Barreled I/A Unit**
- 1410301R **Simcoe Reverse Double-Barreled I/A Unit**



Both models feature .23G. Thin wall irrigation and aspiration. Aspiration through 0.3mm top port with irrigation through side opening, curved shaft, use 28cm of silicone tubing

**1410301** Irrigation through main hub, aspiration through silicone tubing.

**1410301R** Irrigation through silicone tubing, aspiration through main hub.



- 1410302 **Irrigating Vectis** (3x8)



Cannula irrigating vectis, dual fluid flow through to front opening ports

- 1410304 **Vukich Irrigating-Aspirating Cannula**
- 1410304R **Vukich Reverse Irrigating-Aspirating Cannula**



"U" tip irrigation-aspiration cannula designed for manual removal of 12 o'clock codex. 23 G thin wall irrigation and aspiration lines for maximum aspiration capabilities.

**1410304** Irrigation through main hub, aspiration through silicone tubing.

**1410304R** Irrigation through silicone tubing, aspiration through main hub.



ALBARIARTPRINTER  
03004151546

**Posterior Capsule Polisher**

1410303 25 guage  
1410303H 23 guage



Titanium impregnated on posterior surface of olive tip. Front opening. Gently curved

**Air Injection Cannula**

1410305



Angled 45° shaft, 5mm from tip to bend with polished finish, 27G

**Nichamin Hydrodissection Cannula**


1410306



Curved cannula designed to follow the scleral lunel incision and anterior lens capsule to provide less deformation of the wound and less loss of the anterior chamber. The unique flattened tip allows for a better fit under the capsular leaf let. 25G

**Lacrimal Cannula**


1410307



23G. Cannula with blunt tip. Straight shaft.

**Lacrimal Cannula**

1410308



23G. Cannula with reinforced shaft. Straight 10mm tip.

**Lacrimal Cannula**

1410309



23G. Cannula with reinforced shaft. Curved 10mm tip.

**West Lacrimal Cannula**


1410310



23G. Cannula with blunt tip and bottom port. Gently curved shaft

ALBARIARTPRINTER  
0300-4151546

Blade: Sappire  
 Handle: Macromolecular materials



1410051



1.0 x 0.20mm x 15°  
 5.0(mm)  
 Single Edge knife 1.0 x 0.20mm x 15°  
 Length 5.0(mm)

1410053



2.8 x 0.22mm  
 5.0(mm)  
 Tunnel knife 2.8 x 0.22mm  
 Length 5.0(mm)

1410058



3.0 x 0.25mm x 60°  
 7.0(mm)  
 dual-edge kinfe 3.0 x 0.25mm x 60°  
 Length 7.0(mm)

## Sapphire Surgical Knife

Blade: Sappire  
 Handle: Tilanium alloy



1410051



1.0 x 0.20mm x 15°  
 5.0(mm)  
 Single Edge knife 1.0 x 0.20mm x 15°  
 Length 5.0(mm)

1410053



2.8 x 0.22mm  
 5.0(mm)  
 Tunnel knife 2.8 x 0.22mm  
 Length 5.0(mm)

1410058




3.0 x 0.25mm x 60°  
 7.0(mm)  
 dual-edge kinfe 3.0 x 0.25mm x 60°  
 Length 7.0(mm)

ALBARIARTPRINTER  
 0300-4151546

- Moulded high impact resin with see-through lids
- Silicone mats in each tray
- Sterilisation by all standard methods
- Autoclave/Dry Heat 360° F(160°C)/Ethylene Oxide/Cold Solutions


**1410021**



152 x 79 x 17mm (小)

**2~3 Instruments**


**1410022**



254 x 152 x 19mm

**10 Instruments**


ALBARIARTPRINTER  
0300-4151546



- |     |          |   |         |  |  |
|-----|----------|---|---------|--|--|
| 1.  | 1-020S   | Bard Parker Handle                                      | 4-122S* | Halsted Mosquito Forceps, Straight, Long |  |
| 2.  | 1-010T   | Blade Holder  | 4-123S* | Halsted Mosquito Forceps, Curved, Long   |  |
| 3.  | 2-010T   | Castroviejo Caliper                                     | 17.     | 5-060                                    | Gass Retinal Detachment Hook                         |
| 4.  | 4-100S   | Straight Iris Forceps With Serrations                   | 18.     | 5-040                                    | Jameson Muscle Hook                                  |
| 5.  | 4-101S   | Straight Iris Forceps, 1x2 Teeth                        | 19.     | 5-062                                    | Stevens Curved Tenotomy Hook                         |
|     | 4-102S*  | Curved Iris Forceps, 1x2 Teeth                          | 20.     | 6-10/6-050                               | 45 Degree Single Edge Blade                          |
| 6.  | 4-0501T  | Colibri Corneal Forceps, .12 mm 1x2 Teeth, 77 mm        | 21.     | 8-011T*                                  | Barraquer Needle Holder, Small                       |
|     | 4-050T*  | Colibri Corneal Forceps, .12 mm 1x2 Teeth, 84mm         |         | 8-031T                                   | Barraquer Needle Holder, Medium                      |
|     | 4-053T*  | Colibri Corneal Forceps, .12 mm 1x2 Teeth, 115 mm       |         | 8-061T*                                  | Barraquer Needle Holder, Long Size                   |
|     | 4-054T*  | Colibri Corneal Forceps, .12 mm 1x2 Teeth, Round Handle | 22.     | 8-025T*                                  | Barraquer Needle Holder, Extra Fine Jaws, Small      |
|     |          |   |         | 8-045T                                   | Barraquer Needle Holder, Extra Fine Jaws, Medium     |
| 7.  | 4-0600T* | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth          |         | 8-051T*                                  | Barraquer Needle Holder, Extra Fine Jaws, Long       |
|     | 4-058T   | Bonn Corneal Forceps, Straight, .12 mm 1x2 Teeth        | 23.     | 10-021                                   | Desmarres Lid Retractor, Size 1                      |
|     | 4-055T*  | Catalano Corneal Forceps, .12 mm 1x2 Teeth              | 24.     | 10-022*                                  | Desmarres Lid Retractor, Size 2                      |
| 8.  | 4-0601T  | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth           | 25.     | 10-023*                                  | Desmarres Lid Retractor, Size 3                      |
|     | 4-0607S* | Bishop-Harmon Suturing Forceps, .3 mm 1x2 Teeth         | 26.     | 11-040S                                  | Westcott Tenotomy Scissors, Blunt Tips, 13 mm Blades |
| 9.  | 4-0602T* | Castroviejo Suturing Forceps, .5 mm 1x2 Teeth           |         | 11-042S*                                 | Westcott Tenotomy Scissors, Blunt Tips, 6 mm Blades  |
|     | 4-0606T  | Paufique Suturing Forceps, .5 mm 1x2 Teeth              | 27.     | 11-044S                                  | Westcott Stitch Scissors, Sharp Tips, 6 mm Blades    |
| 10. | 4-070S   | Dressing Forceps, 6 mm Serrated Tips                    |         | 11-046S*                                 | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    |
|     | 4-071S*  | Dressing Forceps, 12 mm Serrated Tips                   |         | 11-047S*                                 | Westcott Stitch Scissors, Sharp Tips, 10 mm Blades   |
| 11. | 4-2300T  | Bonaccolto Utility Forceps                              | 28.     | 11-081S                                  | Curved Iris Scissors                                 |
|     | 4-2301T* | Fechtner Conjunctiva Forceps                            | 29.     | 11-080S                                  | Straight Iris Scissors                               |
| 12. | 4-130S   | Jameson Muscle Forceps, Left                            | 30.     | 11-054S                                  | Vannas Capsulotomy Scissors, 6 mm Blades             |
| 13. | 4-131S   | Jameson Muscle Forceps, Right                           |         | 11-0581S*                                | Gills-Vannas Capsulotomy Scissors, 10 mm Blades      |
| 14. | 4-171T*  | McPherson Tying Forceps, 4 mm Platform                  | 31.     | 14-040T                                  | Lieberman Temporal Speculum                          |
|     | 4-178S*  | McPherson Tying Forceps, 7 mm Platform                  | 32.     | 14-022S                                  | Barraquer Wire Speculum, Adult Size                  |
|     | 4-185T   | Tennant Straight Tying Forceps                          | 33.     | 14-023S                                  | Barraquer Wire Speculum, Child Size                  |
|     | 4-186S*  | Tennant Curved Tying Forceps                            | 34.     | 15-051-25                                | Rycroft Anterior Chamber Cannula, 25 Ga              |
| 15. | 4-091S*  | Kelman-McPherson Tying Forceps, 8 mm Platforms          | 35.     | 15-301/303                               | Silicone Bulb With Adapter                           |
|     | 4-092T   | Kelman-McPherson Tying Forceps, 10 mm Platforms         | 36.     | 16-081S                                  | Towel Forceps  |
|     | 4-174T*  | McPherson Angled Tying Forceps, 8 mm Platforms          | 37.     | 16-115                                   | Flynn Scleral Depressor, Overall Length 90 mm        |
|     | 4-175T*  | McPherson Angled Tying Forceps, 10 mm Platforms         | 38.     | 18-304*                                  | Plastic Sterilizing Tray, Single Level, Large        |
| 16. | 4-120S   | Hartman Mosquito Forceps, Straight, Short               | 39.     | 18-305*                                  | Plastic Sterilizing Tray, Double Level, Large        |
|     | 4-121S*  | Hartman Mosquito Forceps, Curved, Short                 |         |  |  |

ALBARIART-PRINTER  
0300-4151546

Extra-Capsular Cataract Extraction (ECCE)




- |     |             |   |     |            |   |
|-----|-------------|---|-----|------------|---|
| 1.  | 2-010T      | Castroviejo Caliper, Titanium                         | 12. | 7-0631*    | Nagahara Phaco Chopper and Drysdale Nucleus Manipulator |
| 2.  | 2-064T      | Adler Wound Gauge, 2.0-2.5 mm                         |     | 7-065      | Rosen Phaco Chopper, Universal                          |
|     | 2-065T*     | Adler Wound Gauge, 2.6-3.0 mm                         | 13. | 7-025      | Ernest Nucleus Cracker                                  |
| 3.  | 3-040       | Whipple Capsulorrhexis Concentration Marker           | 14. | 7-081-23   | Irrigation Handpiece, Curved, 23 Ga                     |
| 4.  | 4-0311T     | Utrata Capsulorrhexis Forceps                         | 15. | 7-0821-23  | Aspiration Handpiece, Curved, 23 Ga                     |
|     | 4-03731     | Kershner One-Pinch Forceps                            | 16. | 8-041T     | Barraquer Needle Holder, Medium                         |
| 5.  | 4-050T*     | Colibri Corneal Forceps, .12 mm 1x2 Teeth             | 17. | 11-040S    | Westcott Curved Tenotomy Scissors                       |
|     | 4-0501T     | Colibri Corneal Forceps, .12 mm 1x2 Teeth             |     | 11-052S*   | Vannas Capsulotomy Scissors                             |
|     | 4-0541T*    | Castroviejo-Colibri Corneal Forceps, .12 mm 1x2 Teeth |     | 11-058S    | Gills-Vannas Capsulotomy Scissors                       |
| 6.  | 4-0600S*    | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth        | 18. | 14-022S    | Barraquer Wire Speculum, Adult Size                     |
|     | 4-0600T     | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth        |     | 14-0221S*  | Barraquer Speculum With Solid Blades, Adult Size        |
|     | 4-0601T*    | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth         |     | 14-025S*   | Barraquer Wire Speculum, Open Blades, Adult Size        |
| 7.  | 4-176T      | McPherson Angled Tying Forceps                        | 19. | 14-040TK   | Lieberman Temporal Speculum                             |
|     | 4-185S*     | Tennant Straight Tying Forceps                        |     | 14-080A*   | Lieberman Temporal Speculum, V Blades, With Aspiration  |
|     | 4-1851T*    | Tennant Angled Tying Forceps                          | 20. | 15-129-0.3 | Simcoe I/A Cannula, 23/23 Ga, 0.3 mm Port, Regular      |
| 8.  | 4-120S      | Hartman Mosquito Forceps, Straight                    |     | 15-0681    | Chang Nucleus Hydrossector                              |
| 9.  | 4-2107T     | Steinert Paddle Lens Folding Forceps                  |     | 15-071-25  | McIntyre Nucleus Hydrossector, Spatulated               |
|     | 4-2141T*    | Cartridge Loading Forceps                             | 21. | 15-301/303 | Silicone Bulb With Adapter                              |
| 10. | 5-030*      | Kuglen Hook, Angled                                   | 22. | 16-080S    | Towel Clamp   |
|     | 5-032       | Sinskey Hook, Angled                                  | 23. | 16-2806    | IOL injector for A, B, C cartridges                     |
|     | 5-0331*     | Lester Lens Manipulator, Angled                       |     | 16-2808*   | IOL injector for D cartridge                            |
|     | 5-034*      | Bechert Nucleus Rotator, Angled                       | 24. | 18-304*    | Plastic Sterilizing Tray, Single Level, Large           |
| 11. | 6-10/6-0501 | 45 Degree Side Port Single Edge Blade 0.6 mm          |     |            |   |
|     | 6-10/6-053* | Trifacet Blade, 1.0 mm                                |     |            |   |
|     | 6-20/6-104* | 2.3/2.8 mm Trapezoid Self-Diving Diamond Knife        |     |            |   |

\* not shown

Set also includes disposable products (see page 199)

## 24-022 CATARACT SET (MINI I)




- |    |         |  |     |            |  |
|----|---------|--|-----|------------|--|
| 1. | 4-050T  | Colibri Corneal Forceps, .12 mm 1x2 Teeth      | 9.  | 8-031T     | Barraquer Needle Holder, Medium Size           |
| 2. | 4-111T  | Jeweler Forceps # 5, Straight                  | 10. | 11-011S    | Castroviejo Universal Corneal Scissors         |
| 3. | 4-136S  | Osher Superior Rectus Forceps                  | 11. | 11-048S    | Westcott Curved Tenotomy Scissors              |
| 4. | 4-2107T | Steinert Paddle Lens Folding Forceps           | 12. | 13-050     | Castroviejo Double Ended Cyclodialysis Spatula |
|    | 4-2108S | Faulkner Lens Holding Forceps                  | 13. | 14-023S    | Barraquer Wire Speculum, Child Size            |
| 5. | 4-0600T | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth | 14. | 15-051-25  | Rycroft Anterior Chamber Cannula, 25 Ga        |
| 6. | 4-0602T | Castroviejo Suturing Forceps, .5 mm 1x2 Teeth  | 15. | 15-129-0.3 | Simcoe I/A Cannula, 23/23 Ga, 0.3 mm Port      |
| 7. | 4-090T  | Kelman-McPherson Tying Forceps                 | 16. | 16-072     | Wilder Lens Loop                               |
| 8. | 5-032   | Sinsky Hook, Angled                            | 17. | 18-304*    | Plastic Sterilizing Tray, Single Level, Large  |

\* not shown

Set also include disposable products (see page 199)


## 24-023 CATARACT SET (MINI II)



- |     |         |  |     |            |   |
|-----|---------|--|-----|------------|---|
| 1.  | 2-010T  | Castroviejo Caliper                            | 12. | 5-040      | Jameson Muscle Hook                             |
| 2.  | 4-050T  | Colibri Corneal Forceps, .12 mm 1x2 Teeth      | 13. | 8-031T     | Barraquer Needle Holder, Medium                 |
| 3.  | 4-136S  | Osher Superior Rectus Forceps                  | 14. | 8-091T     | Barraquer Needle Holder, Long                   |
| 4.  | 4-2107T | Steinert Paddle Lens Folding Forceps           | 15. | 11-011S    | Castroviejo Universal Corneal Scissors          |
|     | 4-2108S | Faulkner Lens Holding Forceps                  | 16. | 11-048S    | Westcott Curved Tenotomy Scissors               |
| 5.  | 4-0600T | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth | 17. | 13-050     | Castroviejo Double Ended Cyclodialysis Spatula  |
| 6.  | 4-0602T | Castroviejo Suturing Forceps, .5 mm 1x2 Teeth  | 18. | 14-023S    | Barraquer Wire Speculum, Child Size             |
| 7.  | 4-111T  | Jeweler Forceps # 5, Straight                  | 19. | 15-051-25  | Rycroft Anterior Chamber Cannula, 25 Ga x 22 mm |
| 8.  | 4-090T  | Kelman-McPherson Tying Forceps                 | 20. | 15-129-0.3 | Simcoe I/A Cannula, 23/23 Ga, 0.3 mm Port       |
| 9.  | 4-121S  | Hartman Mosquito Forceps                       | 21. | 15-301/303 | Silicone Bulb With Adapter                      |
| 10. | 5-032   | Sinsky Hook, Angled                            | 22. | 16-072     | Wilder Lens Loop                                |
| 11. | 5-030   | Kuglen Iris Hook, Angled                       | 23. | 18-304*    | Plastic Sterilizing Tray, Single Level, Large   |

ALBANY PRINTER  
0300-4121546

Phacoemulsification Set




- |     |              |   |     |            |  |
|-----|--------------|---|-----|------------|--|
| 1.  | 2-010T       | Castroviejo Caliper                                   | 21. | 7-0821     | Aspiration Handpiece, Curved, 22Ga, 1 x .35mm Port   |
| 2.  | 4-050T       | Colibri Corneal Forceps, .12 mm 1x2 Teeth             | 22. | 8-031T     | Barraquer Needle Holder, Medium                      |
| 3.  | 4-072T       | Dressing Forceps With 12 mm Serrated Tips             | 23. | 8-045T     | Barraquer Needle Holder, Extra Fine Jaws, Medium     |
| 4.  | 4-2107T      | Steinert Paddle Lens Folding Forceps                  | 24. | 10-5067-1  | Disposable Iris Retractors, 1 pack of 5 pcs          |
|     | 4-2108S      | Faulkner Lens Holding Forceps                         | 25. | 11-040S    | Westcott Tenotomy Scissors, Blunt Tips, 13mm Blades  |
| 5.  | 4-2141T      | Cartridge Loading Forceps                             | 26. | 11-044S    | Westcott Stitch Scissors, Sharp Tips, 16mm Blades    |
| 6.  | 4-0600T      | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth        | 27. | 11-0581S   | Gills-Vannas Capsulotomy Scissors, 10mm Blades       |
|     | 4-0601T*     | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth         | 28. | 11-034S    | Universal Corneal Scissors, Blunt Tips, 7.5mm Blades |
| 7.  | 4-174T       | McPherson Angled Tying Forceps, 8 mm Platforms        | 29. | 14-022S    | Barraquer Wire Speculum, Adult Size                  |
| 8.  | 4-185T       | Tennant Straight Tying Forceps                        | 30. | 14-080A    | Lieberman Temporal Speculum, With Aspiration         |
| 9.  | 4-120S       | Hartman Mosquito Forceps, Straight, Short             | 31. | 14-040T    | Lieberman Temporal Speculum                          |
| 10. | 4-0301T      | Utrata Capsulorhexis Forceps, Cystotome Tips          | 32. | 15-129-0.3 | Simcoe I/A Cannula, 23/23Ga, 0.3mm Port              |
|     | 4-0331T*     | Utrata Capsulorhexis Forceps, Cystotome Tips          | 33. | 15-025     | Kelman Sharp Irrigating Cystotome, 25Ga              |
| 11. | 5-032        | Sinsky Hook, Angled                                   | 34. | 15-079R    | McIntyre-Binkhorst Cannula, 26Ga, Right              |
| 12. | 5-030        | Kuglen Hook, Angled                                   | 35. | 15-079L*   | McIntyre-Binkhorst Cannula, 26Ga, Left               |
| 13. | 5-034        | Bechert Nucleus Rotator, Angled                       | 36. | 15-071-25  | McIntyre Nucleus Hydrossector, Spatulated            |
| 14. | 6-322/6-0531 | 20 Degree Trifacet Blade, 1.0 mm                      | 37. | 15-051-27  | Rycroft Anterior Chamber Cannula, 27Ga               |
| 15. | 6-20/6-091   | 1.50 mm Round Crescent Knife                          | 38. | 15-159-25  | Jensen Capsule Polisher, Curved, 25Ga                |
| 16. | 6-20/6-104*  | 2.3/2.8 mm Trapezoid Self-Diving Diamond Knife        | 39. | 15-301/303 | Silicone Bulb With Adapter                           |
| 17. | 7-070        | Koch Nucleus Spatula                                  | 40. | 16-020T    | Maloney Intraoperative Keratometer                   |
| 18. | 7-063        | Nagahara Phaco Chopper, RHD                           | 41. | 16-081S    | Towel Forceps  |
|     | 7-064*       | Nagahara Phaco Chopper, LHD                           | 42. | 16-0341T   | Fine/Thornton Fixation Ring With Swivel, 14mm        |
| 19. | 7-025T       | Ernest Nucleus Cracker                                | 43. | 16-2806    | IOL injector for A, B, C cartridges                  |
| 20. | 7-081        | Irrigation Handpiece, Curved, 21 Ga, 2 x .35 mm Ports |     | 16-2808*   | IOL injector for D cartridge                         |
|     |              |   | 44. | 18-304*    | Plastic Sterilizing Tray, Single Level, Large        |
|     |              |   | 45. | 18-305*    | Plastic Sterilizing Tray, Double Level, Large        |

ALBARIARTPINTER  
0300-4151546



Phaco Prechop Set




- |     |              |   |     |           |  |
|-----|--------------|---|-----|-----------|--|
| 1.  | 2-010T       | Castroviejo Caliper                                   | 11. | 7-1165S   | Inamura Eagle Prechopper, Straight                   |
| 2.  | 4-0301T      | Utrata Capsulorhexis Forceps, Cystotome Tips          | 12. | 11-040S   | Westcott Tenotomy Scissors, Blunt Tips, 13 mm Blades |
|     | 4-0331T      | Utrata Capsulorhexis Forceps, Cystotome Tips          | 13. | 11-044S   | Westcott Stitch Scissors, Sharp Tips, 16 mm Blades   |
| 3.  | 5-030        | Kuglen Hook, Angled                                   | 14. | 14-040T   | Lieberman Temporal Speculum                          |
| 4.  | 5-032        | Sinskey Hook, Angled                                  | 15. | 15-051-30 | Rycroft Anterior Chamber Cannula, 30 Ga              |
| 5.  | 6-20/6-072   | 2.80 mm Clear Cornea Blade                            | 16. | 15-071-27 | McIntyre Nucleus Hydrossector, 27 Ga                 |
| 6.  | 6-322/6-0531 | 20 Degree Trifacet Blade, 1.0 mm                      | 17. | 15-169-23 | Kratz Capsule Polisher, 23 Ga                        |
| 7.  | 7-025T       | Ernest Nucleus Cracker                                | 18. | 21-R3925  | Irrigating Cystotome - Formed 25 Ga                  |
| 8.  | 7-081        | Irrigation Handpiece, Curved, 21 Ga, 2 x .35 mm Ports | 19. | 16-020T   | Maloney Intraoperative Keratometer                   |
| 9.  | 7-0821       | Aspiration Handpiece, Curved, 22 Ga, 1 x .35 mm Port  | 20. | 18-304*   | Plastic Sterilizing Tray, Single Level, Large        |
| 10. | 7-1161S      | Combo Prechopper, Combo Sharp&Blunt Blades            |     |           |  |

\* not shown

Set also includes disposable products (see page 199)

## Bimanual Microphaco (MICS) Set




1.	4-0501T	Colibri Corneal Forceps, .12 mm 1x2 Teeth	9.	7-063	Nagahara Phaco Chopper, RHD
2.	4-0600T	Castroviejo Suturing Forceps, .12 mm 1x2 Teeth		7-064*	Nagahara Phaco Chopper, LHD
3.	4-174T	McPherson Angled Tying Forceps, 8 mm Platforms	10.	7-070	Koch Nucleus Spatula
4.	4-0374	Capsulorhexis Forceps With 2 Engravings, 23 Ga	11.	7-081	Irrigation Handpiece, Curved, 21 Ga, 2 x .35 mm Ports
	4-0375*	Capsulorhexis Forceps With 6 Engravings, 23 Ga	12.	7-0821	Aspiration Handpiece, Curved, 22 Ga, 1 x .35mm Port
	4-03721	Zaldivar Iridectomy Forceps, With Delicate Tips, 23 Ga	13.	8-045T	Barraquer Needle Holder, Extra Fine Jaws, Medium
	4-03731	Kershner Capsulorhexis Forceps, 23 Ga	14.	11-044S	Westcott Stitch Scissors, Sharp Tips, 16 mm Blades
	4-03741	Kershner Capsulorhexis Forceps, Micro Jaws, 23 Ga	15.	11-03741	Side Port Capsulotomy Scissors, Left, 20 Ga
	4-03751	Fine-Ikeda Capsulorhexis Forceps, Micro Jaws, 23 Ga		11-03751	Side Port Capsulotomy Scissors, Right, 20 Ga
	4-03761*	Ikeda Micro Capsulorhexis Forceps, 23 Ga	16.	14-022S	Barraquer Wire Speculum, Adult Size
	4-03771*	Kawai Capsulorhexis Forceps, 23 Ga	17.	14-080A	Lieberman Temporal Speculum, V Blades, With Aspiration
5.	12-003T*	Instrument Handle (at least 2 pieces)	18.	15-0681	Chang Nucleus Hydrossector
6.	5-032	Sinsky Hook, Angled	19.	15-129-0.3	Simcoe I/A Cannula, 23/23 Ga, 0.3 mm Port
7.	5-034	Bechert Nucleus Rotator, Angled	20.	15-301/303	Silicone Bulb With Adapter
8.	6-20/6-0551	Zaldivar Universal ICL Knife, 0.55/1.00 mm Width	21.	18-304*	Plastic Sterilizing Tray, Single Level, Large
	6-20/6-140*	Trapezoid Side Port Knife, 1.00 mm	22.	18-305*	Plastic Sterilizing Tray, Double Level, Large
	6-20/6-142	Micro Incision Knife, 1.30/1.50 mm			
	6-20/6-144*	Micro Incision Knife, 1.80/2.00 mm			

\* not shown

Set also includes disposable products (see page 199)

LASIK Set



- |    |            |   |          |  |  |
|----|------------|---|----------|--|--|
| 1. | 3-174T     | Lavery LASIK Marker                             | 15-378R  | Gimbel LASIK Fountain Cannula, Right, 25Ga |  |
| 2. | 3-176T     | LASIK Flap Marker                               | 15-378L* | Gimbel LASIK Fountain Cannula, Left, 25Ga  |  |
| 3. | 13-110     | Paton Double Ended Spatula And Spoon            | 6.       | 15-376                                     | Slade LASIK Cannula, 26G                         |
| 4. | 14-040TL   | Lieberman Temporal Speculum, Rounded Blades     | 7.       | 15-379-25                                  | Buratto LASIK Irrigation Cannula 25Ga            |
|    | 14-0401TL* | Lieberman Temporal Speculum, Flatten Blades     | 8.       | 16-0341T                                   | Fine/Thornton Fixation Ring With Swivel, 14mm    |
|    | 14-0411TL* | Lieberman Nasal Speculum, Flatten Blades, Short | 9.       | 18-304*                                    | Plastic SterilizingTray, Single Level, Large     |
|    | 14-080LA   | Lieberman Temporal Speculum                     | 10.      | 20-001                                     | Hockey Epithelium Removal Knife                  |
|    | 14-081LA*  | Lieberman Nasal Speculum                        | 11.      | 20-002                                     | Lindstrom LASIK Spatula&Epithelium Removal Board |
| 5. | 15-373-25  | Banaji LASIK Irrigation Cannula, 25Ga           | 12.      | 20-013                                     | LASIK Spatula And Flap Retreatment Instrument    |
|    | 15-373-27* | Banaji LASIK Irrigation Cannula, 27Ga           | 13.      | 20-050                                     | Disposable Marking Pen, pack of 10               |

\* not shown

24-032 REFRACTIVE SET

LASEK Set




- |    |            |   |     |           |   |
|----|------------|---|-----|-----------|---|
| 1. | 14-040TL   | Lieberman Temporal Speculum, Rounded Blades     | 5.  | 20-1011   | LASEK Trephine 8.0mm, Preset Depth 70 Microns |
|    | 14-0401TL* | Lieberman Temporal Speculum, Flatten Blades     |     | 20-1021*  | LASEK Trephine 9.0mm, Preset Depth 70 Microns |
|    | 14-0411TL* | Lieberman Nasal Speculum, Flatten Blades, Short | 6.  | 20-1031T  | LASEK Funnel 8.5 mm                           |
|    | 14-080LA   | Lieberman Temporal Speculum                     |     | 20-1041T* | LASEK Funnel 9.5 mm                           |
|    | 14-081LA*  | Lieberman Nasal Speculum                        | 7.  | 20-130    | LASEK Detaching Spatula                       |
| 2. | 15-097     | Troutman Alpha-Chrymostrypsin Cannula, 25Ga     | 8.  | 20-131    | LASEK Epithelial Micro Hoe                    |
| 3. | 20-121     | LASEK 8.0mm Trephine / 8.5mm Funnel             | 9.  | 20-0011   | LASEK Knife                                   |
| 4. | 20-122     | LASEK 9mm Trephine / 9.5mm Funnel               | 10. | 20-110    | LASEK Spatula                                 |
|    |            |   | 11. | 18-304    | Plastic SterilizingTray, Single Level, Large  |

\* not shown

ALBARIARTPRINTER  
0300-4151546

FemtoLASIK Set




- |    |            |   |    |           |  |
|----|------------|---|----|-----------|--|
| 1. | 3-174T     | Lavery LASIK Marker                             | 3. | 14-080LA  | Lieberman Temporal Speculum                  |
|    | 3-176T     | LASIK Flap Marker                               |    | 14-081LA* | Lieberman Nasal Speculum                     |
| 2. | 14-040TL   | Lieberman Temporal Speculum, Rounded Blades     | 4. | 20-201    | FemtoLASIK Flap Spatula Double Ended         |
|    | 14-0401TL* | Lieberman Temporal Speculum, Flatten Blades     | 5. | 20-202    | FemtoLASIK Flap Spatula                      |
|    | 14-0411TL  | Lieberman Nasal Speculum, Flatten Blades, Short | 6. | 20-203    | Zaldivar FemtoLASIK Spatula                  |
|    |            |   | 7. | 15-379-25 | Buratto LASIK Irrigation Cannula 25Ga        |
|    |            |   | 8. | 18-304    | Plastic SterilizingTray, Single Level, Large |

\* not shown

24-034 REFRACTIVE SET

ICL™ Set




- |    |             |   |     |         |  |
|----|-------------|---|-----|---------|--|
| 1. | 4-20111S    | ICL™ Cartridge Loading Forceps                    | 6.  | 12-003T | Instrument Handle                            |
| 2. | 4-21431     | Coaxial Angled ICL™ Loading Forceps, 20 Ga        | 7.  | 13-141  | Pallikaris ICL™ Manipulator                  |
| 3. | 4-08011T    | Nevyas-Wallace Fixation Forceps                   | 8.  | 13-142  | Zaldivar ICL™ Manipulator                    |
| 4. | 6-20/6-104  | 2.3/2.8mm Trapezoid Self-Diving Diamond Knife     | 9.  | 14-040T | Lieberman Temporal Speculum                  |
| 5. | 6-20/6-0551 | Zaldivar Universal ICL™ Knife, 0.55/1.00 mm Width | 10. | 18-304  | Plastic SterilizingTray, Single Level, Large |

ICL™ is a trade mark of STAAR Surgical Company

ALBARIARTPRINTER  
0300-4151546

LRI Set




- |    |              |  |                |   |  |
|----|--------------|--|----------------|---|--|
| 1. | 2-0311T      | LRI Gauge with atraumatic fixation teeth                                       | 6-500/6-0531   | 20 Degree Trifacet Blade, 0,20/1.0 mm<br>0.50 mm preset depth               |  |
|    | 2-0331T      | Mendez Degree Gauge with grooves   | 6-500S/6-0531  | 21 Degree Trifacet Blade, 0,20/1.0 mm<br>0.50 mm preset depth, short handle |  |
| 2. | 3-180        | LRI Marker, 6.0&8.0mm Chord Lengths  | 6-600/6-0531*  | 20 Degree Trifacet Blade, 0,20/1.0 mm<br>0.60 mm preset depth               |  |
|    | 3-1801       | LRI Marker, 40-60-80 degrees   | 6-600S/6-0531* | 21 Degree Trifacet Blade, 0,20/1.0 mm<br>0.60 mm preset depth, shdrandle    |  |
| 3. | 3-193        | Whitehouse Gravity Axis Marker, Angled   | 6.             | 14-040T   | Lieberman Temporal Speculum                  |
|    | 3-1931*      | Whitehouse Gravity Axis Marker, Straigh  | 7.             | 18-304  | Plastic SterilizingTray, Single Level, Large |
|    | 3-1932       | Whitehouse Gravity Axis Marker II, Straight,<br>Reduced                        | 8.             | 20-050  | Disposable Marking Pen, pack of 10           |
| 4. | 4-08011T     | Nevyas-Wallace Fixation Forceps  |                |   |  |
| 5. | 6-322/6-0531 | 20 Degree Trifacet Blade, 1.0 mm width, 0.50,<br>0.55 & 0.60 mm incision depth |                |   |  |

\* not shown

24-036 REFRACTIVE SET

TORIC IOL Set



- |    |         |   |    |         |  |
|----|---------|---|----|---------|--|
| 1. | 2-0311T | LRI Gauge with atraumatic<br>fixation teeth             | 3. | 3-194   | Rumex Toric Combo Marker, Vertical 0-0 Axis      |
|    | 2-0331T | Mendez Degree Gauge with grooves                        |    | 3-1941  | Rumex Toric Combo Marker II, Horizontal 0-0 Axis |
| 2. | 3-193   | Whitehouse Gravity Axis Marker, Angled                  | 4. | 3-090T  | Bores Axis Marker                                |
|    | 3-1931* | Whitehouse Gravity Axis Marker, Straigh                 | 5. | 14-040T | Lieberman Temporal Speculum                      |
|    | 3-1932  | Whitehouse Gravity Axis Marker II,<br>Straight, Reduced | 6. | 15-170  | MICS Capsule Polisher, 23 Ga                     |
|    |         |   | 7. | 18-304  | Plastic SterilizingTray, Single Level, Large     |
|    |         |   | 8. | 20-050  | Disposable Marking Pen, pack of 10               |

\* not shown

ALBARIART-PRINTER  
0300-4151546


## CORNEAL TRANSPLANTATION Set



- |     |            |   |     |            |   |
|-----|------------|---|-----|------------|---|
| 1.  | 3-140T     | Corneal Transplant Marker                             | 15. | 11-024S    | Castroviejo Corneal Section Scissors, Right   |
| 2.  | 4-0814T    | Pollack Corneal Transplantation Fixation Forceps      | 16. | 11-0241S*  | Castroviejo Corneal Section Scissors, Left    |
| 3.  | 4-0541T    | Castroviejo Colibri Corneal Forceps, .12 mm 1x2 Teeth | 17. | 11-040S    | Westcott Curved Tenotomy Scissors             |
| 4.  | 4-0607S    | Bishop-Harmon Suturing Forceps, .3 mm 1x2 Teeth       | 18. | 11-0581S   | Gills-Vannas Capsulotomy Scissors             |
| 5.  | 4-0600T    | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth        | 19. | 11-044S*   | Westcott Stitch Scissors                      |
| 6.  | 4-0601T*   | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth         | 20. | 13-110     | Paton Double Ended Spatula And Spoon          |
| 7.  | 4-090T     | Kelman-McPherson Tying Forceps, 4 mm Platforms        | 21. | 14-022S    | Barraquer Wire Speculum, Adult Size           |
| 8.  | 4-178S     | McPherson Straight Tying Forceps, 7mm Platforms       | 22. | 15-051-25  | Rycroft Anterior Chamber Cannula, 25 Ga       |
| 9.  | 4-120S     | Hartman Mosquito Forceps                              | 23. | 15-301/303 | Silicone Bulb With Adapter                    |
| 10. | 6-10/6-056 | Alfonso Corneal Transplant Knife                      | 24. | 16-081S*   | Towel Forceps                                 |
| 11. | 8-031T     | Barraquer Needle Holder, Medium                       | 25. | 16-020T    | Maloney Intraoperative Keratometer            |
| 12. | 8-045T     | Barraquer Needle Holder, Extra Fine Jaws, Medium      | 26. | 16-030S*   | Corneal Trephine Blades, 7.50 mm              |
| 13. | 11-020S    | Katzin Corneal Transplant Scissors, Left              | 27. | 16-0307    | Corneal Trephine Blades, 8.0 mm               |
| 14. | 11-0201S*  | Katzin Corneal Transplant Scissors, Right             | 28. | 18-305*    | Plastic Sterilizing Tray, Double Level, Large |


## 24-042 CORNEAL SETS

## DSEK, DSAEK, DMEK Set



- |    |          |  |         |                                     |  |
|----|----------|--|---------|-------------------------------------|--|
| 1. | 3-0219T  | Hoffer Optical Zone Marker, 9.00mm           | 9.      | 13-150T                             | Spatula-Guide for Corneal Endothelium Implantation |
| 2. | 4-034    | Forceps for Corneal Endothelium Implantation | 10.     | 13-151S                             | Cindy Sweeper DSEK Spatula                         |
| 3. | 4-2019T  | Corneal Donor Insertion Forceps              | 11.     | 13-152S                             | Carlson DSEK Smoother                              |
| 4. | 5-0322   | Reversed Sinsky Hook                         | 12.     | 13-153S                             | Terry DSEK Scraper                                 |
| 5. | 12-003T  | Instrument Handle                            | 13-154T | Melles DSAEK PLK Scraper, 45 Degree |  |
| 6. | 13-137   | Corneal Dissector, Straight                  | 13-155T | Melles DSAEK PLK Scraper, 90 Degree |  |
| 7. | 13-138   | Corneal Dissector, Curved                    | 13.     | 18-304                              | Plastic Sterilizing Tray, Single Level, Large      |
| 8. | 13-139/I | Endothelial Stripper, Irrigating             |         |                                     |  |


DALK Set



- |    |         |  |    |           |   |
|----|---------|--|----|-----------|---|
| 1. | 13-170  | Trisector for DALK Procedure             | 5. | 11-02015  | Katzin Corneal Transplant Scissors, Right     |
| 2. | 13-171  | Spatula for DALK Procedure               | 6. | 15-450-27 | Cannula for DALK Procedure, 27 Ga             |
| 3. | 13-172  | Dissector for DALK Procedure             | 7. | 18-304    | Plastic Sterilizing Tray, Single Level, Large |
| 4. | 11-0205 | Katzin Corneal Transplant Scissors, Left |    |           |   |

24-044 CORNEAL SETS

DLEK Set



- |    |          |  |     |         |   |
|----|----------|--|-----|---------|---|
| 1. | 2-010T   | Castroviejo Caliper                              | 6.  | 11-0365 | DLEK Scissors, Medium Curve                   |
| 2. | 3-0304T  | Osher-Neuman Low Profile Corneal Marker, 8 Lines | 7.  | 13-137  | Corneal Dissector, Straight                   |
| 3. | 4-0814T  | Pollack Corneal Transplantation Fixation Forceps | 8.  | 13-138  | Corneal Dissector, Curved                     |
| 4. | 5-0322   | Reversed Sinskey Hook                            | 9.  | 13-160  | Manipulator for DLEK procedure, Z-hook        |
| 5. | 11-0361S | DLEK Scissors, Strong Curve                      | 10. | 18-304  | Plastic Sterilizing Tray, Single Level, Large |

ALBARIAS PRINTING  
0300-4151546

## 24-05 CHALAZION SET




- |    |        |   |     |         |   |
|----|--------|---|-----|---------|---|
| 1. | 1-020S | Bard Parker Handle                            | 7.  | 4-1906T | Desmarres Chalazion Forceps, Large            |
| 2. | 16-063 | Meyerhoefer Chalazion Curette, Size 0-1.5 mm  | 8.  | 4-1912T | Desmarres Chalazion Forceps, Medium           |
| 3. | 16-064 | Meyerhoefer Chalazion Curette, Size 1-1.75 mm | 9.  | 4-1908T | Lambert Chalazion Forceps, Small              |
| 4. | 16-065 | Meyerhoefer Chalazion Curette, Size 2-2.0 mm  | 10. | 4-1909T | Lambert Chalazion Forceps, Medium             |
| 5. | 16-066 | Meyerhoefer Chalazion Curette, Size 3-2.5 mm  | 11. | 11-044S | Westcott Stitch Scissors, Sharp Tips          |
| 6. | 16-067 | Meyerhoefer Chalazion Curette, Size 4-3.5 mm  | 12. | 18-304  | Plastic Sterilizing Tray, Single Level, Large |

## 24-06 ENUCLEATION/EVISCERATION SET




- |    |         |  |     |         |   |
|----|---------|--|-----|---------|---|
| 1. | 4-0741S | Adson Fixation Forceps, .5 mm 1x2 Teeth        | 7.  | 11-100S | Knapp Straight Strabismus Scissors            |
| 2. | 4-120S  | Hartman Mosquito Forceps, Straight, Short      |     | 11-101S | Knapp Curved Strabismus Scissors              |
| 3. | 4-0540T | Colibri-Bonn Corneal Forceps, .12 mm 1x2 Teeth | 8.  | 14-040T | Lieberman Temporal Speculum                   |
| 4. | 5-042   | Graefe Muscle Hook, Size 2                     | 9.  | 16-060  | Wells Enucleation Spoon                       |
| 5. | 11-090S | Curved Enucleation Scissors, Short             | 10. | 16-061  | Bunge Evisceration Spoon, Small               |
| 6. | 11-091S | Curved Enucleation Scissors, Long              | 11. | 16-062  | Bunge Evisceration Spoon, Large               |
|    |         |  | 12. | 18-304* | Plastic Sterilizing Tray, Single Level, Large |






- |     |            |   |     |            |   |
|-----|------------|---|-----|------------|---|
| 1.  | 1-010T     | Blade Holder                                  | 15. | 11-062S    | McPherson-Vannas Curved Iris Scissors                   |
| 2.  | 2-010T     | Castroviejo Caliper                           | 16. | 11-080S    | Iris Scissors, Straight                                 |
| 3.  | 4-050T     | Colibri Corneal Forceps, .12mm, 1x2 Teeth     | 17. | 14-022S    | Barraquer Wire Speculum, Adult Size                     |
| 4.  | 4-070S     | Dressing Forceps With Delicate Serrations     | 18. | 16-010     | Rumex Corneoscleral Punch<br>(.5, .75, 1.0, 1.5mm Tips) |
| 5.  | 4-0600T    | Castroviejo Suturing Forceps, .12mm 1x2 Teeth |     | 16-011     | Kelly Descemet's Membrane Punch, .75mm Diameter         |
| 6.  | 4-171S     | McPherson Straight Tying Forceps              | 19. | 16-012S    | Harms Trabeculotome, Left                               |
| 7.  | 4-090T     | Kelman-McPherson Tying Forceps                | 20. | 16-013S    | Harms Trabeculotome, Right                              |
| 8.  | 4-120S     | Hartman Mosquito Forceps, Straight            | 21. | 16-081S    | Towel Forceps   |
| 9.  | 6-10/6-050 | 45 Degree Single Edge Blade/Straight Handle   | 22. | 15-051-27  | Rycroft Anterior Chamber Cannula, 27Ga                  |
| 10. | 6-20/6-092 | 2.0mm Round Crescent Knife/Angled Handle      | 23. | 15-301/303 | Silicone Bulb With Adapter                              |
| 11. | 8-031T     | Barraquer Needle Holder, Standard Jaws        | 24. | 18-305*    | Plastic Sterilizing Tray, Double Level, Large           |
| 12. | 8-045T     | Barraquer Needle Holder, Extra Fine Jaws      | 25. | 18-301-1*  | Plastic Sterilizing Tray For 2 Diamond Knives           |
| 13. | 11-040S    | Westcott Curved Tenotomy Scissors             | 26. | 20-050     | Disposable Marking Pen, pack of 10                      |
| 14. | 11-044S    | Westcott Stitch Scissors                      |     |            |   |



- |    |            |                                 |     |         |   |
|----|------------|---------------------------------|-----|---------|---|
| 1. | 2-0311T    | Mendez Degree Gauge             | 7.  | 10-035  | Elevator for ICSR implantation                |
| 2. | 3-034      | Optical Zone Marker             | 8.  | 13-146  | Suarez Spreader                               |
| 3. | 3-143T     | Tunnel Zone Marker              | 9.  | 13-147  | Bicalto Guide                                 |
| 4. | 4-08011T   | Neyyas-Wallace Fixation Forceps | 10. | 16-173S | Tunnel Maker Left                             |
| 5. | 4-2144T    | Forceps for ICSR Implant        | 11. | 16-174S | Tunnel Maker Right                            |
| 6. | 6-00/6-020 | Diamond Knife                   | 12. | 18-304  | Plastic Sterilizing Tray, Single Level, Large |



- |     |          |  |     |             |   |
|-----|----------|--|-----|-------------|---|
| 1.  | 1-010T   | Blade Holder                                     | 13. | 9-060T      | Castroviejo Double-Ended Lacrimal Dilator, Size 1 & 2 |
| 2.  | 4-071S   | Dressing Forceps, 12 mm Serrated Tips            | 14. | 10-013      | Stevenson Lacrimal Sac Retractor                      |
| 3.  | 4-2300T  | Bonaccolto Utility Forceps                       | 15. | 10-014      | Knapp Lacrimal Sac Retractor                          |
| 4.  | 4-0601T  | Castroviejo Suturing Forceps, 3 mm 1x2 Teeth     | 16. | 11-080S     | Straight Iris Scissors                                |
| 5.  | 4-08011T | Neyvas-Wallace Fixation Forceps                  |     | 11-081S*    | Curved Iris Scissors                                  |
| 6.  | 4-120S   | Hartman Mosquito Forceps, Straight, Short        | 17. | 11-125S     | Westcott Type Stitch Scissors                         |
|     | 4-121S*  | Hartman Mosquito Forceps, Curved, Short          | 18. | 15-032      | Shahinian Lacrimal Cannula, 25 Ga                     |
| 7.  | 8-031T   | Barraquer Needle Holder, Medium                  | 19. | 15-029      | Lacrimal Cannula, Reinforced, Curved                  |
| 8.  | 8-045T   | Barraquer Needle Holder, Extra Fine Jaws, Medium |     | 15-031      | Lacrimal Cannula (Bailey)                             |
| 9.  | 9-010S   | Bowman Lacrimal Probe, Size 0000-000             | 20. | 15-033      | Lacrimal Cannula (Fasanella)                          |
|     | 9-011S   | Bowman Lacrimal Probe, Size 00-0                 |     | 15-035-25C* | Lacrimal Cannula (Anel) Curved                        |
|     | 9-012S*  | Bowman Lacrimal Probe, Size 1-2                  | 21. | 15-051-23   | Rycroft Anterior Chamber Cannula, 23 Ga               |
|     | 9-013S*  | Bowman Lacrimal Probe, Size 3-4                  | 22. | 15-301/303  | Silicone Bulb With Adapter                            |
|     | 9-014S*  | Bowman Lacrimal Probe, Size 5-6                  | 23. | 16-081S     | Towel Forceps   |
|     | 9-015S*  | Bowman Lacrimal Probe, Size 7-8                  | 24. | 16-135      | Surgical Mallet, Polished Finish                      |
| 10. | 9-021S   | Quickert Lacrimal Intubation Probe, Size 0       | 25. | 16-136      | Kerrison Ronger, Size 30, mm Wide                     |
| 11. | 9-031    | Pigtail Lacrimal Probe                           | 26. | 18-305*     | Plastic Sterilizing Tray, Double Level, Large         |
| 12. | 9-050T*  | Wilder Lacrimal Dilator, Size 1                  | 27. | 21-R7011    | DCR Set Straight 4.5 cm                               |
|     | 9-051T*  | Wilder Lacrimal Dilator, Size 2                  |     | 21-R7012*   | DCR Set Angled 4.5 cm                                 |
|     | 9-052T   | Wilder Lacrimal Dilator, Size 3                  |     |             |   |


\* not shown

ALBARIARTPRINTER  
0300-4151546




- |     |          |  |     |            |   |
|-----|----------|--|-----|------------|---|
| 1.  | 1-010T   | Blade Holder                                   | 16. | 10-020     | Desmarres Lid Retractor, Size 0                     |
| 2.  | 2-010T   | Castroviejo Caliper                            | 17. | 10-021     | Desmarres Lid Retractor, Size 1                     |
| 3.  | 4-0600T  | Castroviejo Suturing Forceps, .12mm 1x2 Teeth  | 18. | 10-022*    | Desmarres Lid Retractor, Size 2                     |
|     | 4-0601T* | Castroviejo Suturing Forceps, .3mm 1x2 Teeth   | 19. | 10-023*    | Desmarres Lid Retractor, Size 3                     |
| 4.  | 4-0607S  | Bishop-Harmon Suturing Forceps, .3mm 1x2 Teeth | 20. | 11-040S    | Westcott Tenotomy Scissors, Blunt Tips, 13mm Blades |
| 5.  | 4-072T   | Dressing Forceps With 12mm Serrated Tips       | 21. | 11-044S    | Westcott Stitch Scissors, Sharp Tips, 16mm Blades   |
| 6.  | 4-090T   | Kelman-McPherson Tying Forceps                 |     | 11-046S*   | Westcott Stitch Scissors, Sharp Tips, 13mm Blades   |
| 7.  | 4-120S   | Hartman Mosquito Forceps, Straight, Short      | 22. | 11-100S    | Knapp Straight Strabismus Scissors                  |
| 8.  | 4-121S   | Hartman Mosquito Forceps, Curved, Short        |     | 11-101S    | Knapp Curved Strabismus Scissors                    |
| 9.  | 4-140T   | Putterman Type Lid Clamp                       | 23. | 14-040T    | Lieberman Temporal Speculum                         |
| 10. | 4-1906T  | Desmarres Chalazion Forceps, Large             | 24. | 15-301/303 | Silicone Bulb With Adapter                          |
|     | 4-1908T* | Lambert Chalazion Forceps, Small               | 25. | 16-081S    | Towel Forceps                                       |
|     | 4-1912T* | Desmarres Chalazion Forceps, Medium            | 26. | 16-50S     | Lid Plate, Stainless Steel                          |
| 11. | 4-2300T  | Bonaccolto Utility Forceps                     |     | 16-50T     | Lid Plate, Titanium                                 |
| 12. | 5-042    | Graefe Muscle Hook, Size 2                     | 27. | 18-304*    | Plastic Sterilizing Tray, Single Level, Large       |
| 13. | 8-031T   | Barraquer Needle Holder, Medium                | 28. | 18-305*    | Plastic Sterilizing Tray, Double Level, Large       |
| 14. | 8-080T   | Kalt Needle Holder                             | 29. | 20-050     | Disposable Marking Pen, pack of 10                  |
| 15. | 10-014   | Knapp Lacrimal Sac Retractor                   |     |            |   |

\* not shown



- |     |          |   |     |            |   |
|-----|----------|---|-----|------------|---|
| 1.  | 1-020S   | Bard Parker Handle                            | 15. | 5-062      | Stevens Curved Tenotomy Hook                        |
| 2.  | 2-010T   | Castroviejo Caliper                           | 16. | 8-031T     | Barraquer Needle Holder, Medium                     |
| 3.  | 4-072T   | Dressing Forceps With 12mm Serrated Tips      | 17. | 8-091T     | Barraquer Needle Holder, Long                       |
| 4.  | 4-0600T  | Castroviejo Suturing Forceps, .12mm 1x2 Teeth | 18. | 10-020     | Desmarres Lid Retractor, Size 0                     |
|     | 4-0602T* | Castroviejo Suturing Forceps, .5mm 1x2 Teeth  |     | 10-021*    | Desmarres Lid Retractor, Size 1                     |
| 5.  | 4-08011T | Nevyas-Wallace Fixation Forceps               |     | 10-022*    | Desmarres Lid Retractor, Size 2                     |
| 6.  | 4-090T   | Kelman-McPherson Tying Forceps                | 19. | 11-040S    | Westcott Tenotomy Scissors, Blunt Tips, 13mm Blades |
| 7.  | 4-120S   | Hartman Mosquito Forceps, Straight, Short     | 20. | 11-044S    | Westcott Stitch Scissors, Sharp Tips, 16mm Blades   |
|     | 4-121S   | Hartman Mosquito Forceps, Curved, Short       | 21. | 11-080S    | Straight Iris Scissors                              |
| 8.  | 4-130S   | Jameson Muscle Forceps, Left                  |     | 11-081S*   | Curved Iris Scissors                                |
|     | 4-131S*  | Jameson Muscle Forceps, Right                 | 22. | 14-0601T   | Kershner Reversible Solid Blade Speculum            |
| 9.  | 4-136S   | Osher Superior Rectus Forceps                 | 23. | 15-051-23  | Rycroft Anterior Chamber Cannula, 23 Ga             |
| 10. | 5-040    | Jameson Muscle Hook, 2.0mm Bulbous Tip        |     | 15-051-25* | Rycroft Anterior Chamber Cannula, 25 Ga             |
|     | 5-0401*  | Jameson Muscle Hook, 1.5mm Bulbous Tip        |     | 15-051-27* | Rycroft Anterior Chamber Cannula, 27 Ga             |
| 11. | 5-041*   | Graefe Muscle Hook, Size 1                    |     | 15-051-30* | Rycroft Anterior Chamber Cannula, 30 Ga             |
| 12. | 5-042    | Graefe Muscle Hook, Size 2                    | 24. | 15-301/303 | Silicone Bulb With Adapter                          |
| 13. | 5-060    | Gass Retinal Detachment Hook                  | 25. | 16-081S    | Towel Forceps                                       |
| 14. | 5-061*   | Scobee Oblique Muscle Hook, Curved            | 26. | 18-305*    | Plastic Sterilizing Tray, Double Level, Large       |

\* not shown



- |     |          |   |     |          |   |
|-----|----------|---|-----|----------|---|
| 1.  | 1-020T   | Blade Holder                              | 11. | 9-050T   | Wilder Lacrimal Dilator, Size 1               |
|     | 1-020S*  | Bard Parker Handle                        |     | 9-051S*  | Wilder Lacrimal Dilator, Size 2               |
| 2.  | 4-120S   | Hartman Mosquito Forceps, Straight, Short |     | 9-051T*  | Wilder Lacrimal Dilator, Size 2               |
|     | 4-121S*  | Hartman Mosquito Forceps, Curved, Short   |     | 9-052S*  | Wilder Lacrimal Dilator, Size 3               |
| 3.  | 4-140T   | Putterman Type Lid Clamp                  |     | 9-052T*  | Wilder Lacrimal Dilator, Size 3               |
| 4.  | 4-1906T  | Desmarres Chalazion Forceps, Large        | 12. | 10-014   | Knapp Lacrimal Sac Retractor                  |
|     | 4-1908T* | Lambert Chalazion Forceps, Small          | 13. | 10-022   | Desmarres Lid Retractor, Size 2               |
|     | 4-1912T* | Desmarres Chalazion Forceps, Medium       | 14. | 11-090S  | Curved Enucleation Scissors, Short            |
| 5.  | 4-0741S  | Adson Fixation Forceps, .5mm 1x2 Teeth    |     | 11-091S* | Curved Enucleation Scissors, Long             |
| 6.  | 5-042    | Graefe Muscle Hook, Size 2                | 15. | 11-100S  | Knapp Straight Strabismus Scissors            |
| 7.  | 8-031T   | Barraquer Needle Holder, Medium           |     | 11-101S* | Knapp Curved Strabismus Scissors              |
| 8.  | 8-080T   | Kalt Needle Holder                        | 16. | 14-040T  | Lieberman Temporal Speculum                   |
| 9.  | 9-012S   | Bowman Lacrimal Probe, Size 1-2           | 17. | 16-50S   | Lid Plate, Stainless Steel                    |
|     | 9-013S*  | Bowman Lacrimal Probe, Size 3-4           |     | 16-50T*  | Lid Plate, Titanium                           |
|     | 9-014S*  | Bowman Lacrimal Probe, Size 5-6           | 18. | 16-061   | Bunge Evisceration Spoon, Small               |
|     | 9-015S*  | Bowman Lacrimal Probe, Size 7-8           |     | 16-062*  | Bunge Evisceration Spoon, Large               |
| 10. | 9-031    | Pigtail Lacrimal Probe                    | 19. | 16-066   | Meyerhoefer Chalazion Curette, Size 3-2.5mm   |
|     |          |   | 20. | 18-304*  | Plastic Sterilizing Tray, Single Level, Large |
|     |          |   | 21. | 18-305*  | Plastic Sterilizing Tray, Double Level, Large |

\* not shown


## 24-13 PTERYGIUM SET



- |    |         |  |     |             |  |
|----|---------|--|-----|-------------|--|
| 1. | 1-010T  | Blade Holder                                 | 8.  | 14-022S     | Barraquer Wire Speculum, Adult Size                |
| 2. | 4-0602T | Castroviejo Suturing Forceps, .5mm 1x2 Teeth | 9.  | 16-051-2.5  | Power Handle with 2.5 mm Round Diamond Dusted Burr |
| 3. | 4-171T  | McPherson Straight Tying Forceps             | 10. | 16-051-3.5B | Round Diamond Dusted Burr 3.5 mm Diameter          |
| 4. | 4-2300T | Bonaccolto Conjunctiva Forceps               | 11. | 16-052-5.0B | Disk Shaped Diamond Dusted Burr 5.0 mm Diameter    |
| 5. | 4-0814T | Pollack Fixation Forceps                     | 12. | 18-304*     | Plastic Sterilizing Tray, Single Level, Large      |
| 6. | 8-031T  | Barraquer Needle Holder, Curved              |     |             |  |
| 7. | 11-080S | Straight Iris Scissors                       |     |             |  |

\* not shown

ALBARRAN PRINTING  
0300-4151546



- |     |            |  |     |                 |  |
|-----|------------|--|-----|-----------------|--|
| 1.  | 1-020S     | Bard Parker Handle                                   | 23. | 16-111S         | Schocket Double End Cannula                                    |
| 2.  | 1-010T     | Blade Holder   | 24. | 18-304*         | Plastic Sterilization Tray, Single Level, Large                |
| 3.  | 2-010T     | Castroviejo Caliper                                  | 25. | 18-305*         | Plastic Sterilization Tray, Double Level, Large                |
| 4.  | 2-100S     | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 26. | 18-307*         | Mini Sterilization Tray with Silicone Mat, 38 x 68,6 x 25,4 mm |
|     | 2-100T*    | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 27. | 12-001T**       | Vitreoretinal Instrument Handle, One Finger Control            |
|     | 2-101T*    | Braunstein Fixed Caliper, 3.0-3.5 mm                 |     | 12-003T**       | Vitreoretinal Instrument Handle, Squeeze Model                 |
| 5.  | 4-0600T    | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth       | 28. | 12-202          | Vertical Vitreoretinal Scissors, 70 Degree                     |
| 6.  | 4-0601T    | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth        |     | 12-2020*        | Claus Lucke Retinotomy Scissors with Bulbous Tip               |
| 7.  | 4-0607S    | Bishop-Harmon Suturing Forceps, 3 mm 1x2 Teeth       | 29. | 12-209          | Curved Subretinal Scissors, Curvature Radius 12 mm             |
| 8.  | 4-091S     | Kelman-McPherson Tying Forceps, 8 mm Platforms       | 30. | 12-301          | Diamondized End-Gripping Vitreoretinal Forceps                 |
|     | 4-092T*    | Kelman-McPherson Tying Forceps, 10 mm Platforms      |     | 12-304*         | Crocodile Vitreoretinal Gripping Forceps                       |
|     | 4-174T*    | McPherson Angled Tying Forceps, 8 mm Platforms       |     | 12-3044*        | Lucke Multipurpose Forceps                                     |
|     | 4-175T*    | McPherson Angled Tying Forceps, 10 mm Platforms      |     | 12-401*         | Vitreoretinal End Gripping Forceps With Micro Jaws             |
| 9.  | 4-2201T    | Watzke Sleeve Spreading Forceps                      |     | 12-402*         | Vitreoretinal ILM «Nail» End Gripping Forceps                  |
| 10. | 4-2300T    | Bonaccolto Utility Forceps                           | 31. | 12-410          | Eckardt End-Gripping Forceps                                   |
| 11. | 4-122S     | Halsted Mosquito Forceps, Straight, Long             |     | 12-411*         | Tano Asymmetrical End-Gripping Forceps                         |
|     | 4-123S*    | Halsted Mosquito Forceps, Curved, Long               | 32. | 12-335          | Stolyarenko Forceps For Large Foreign Body                     |
| 12. | 5-042      | Graefe Muscle Hook, Size 2                           |     | 12-321*         | Spring Gripping Forceps  |
| 13. | 5-060      | Gass Retinal Detachment Hook                         |     | 12-313*         | Vitreoretinal Forceps With Cup Jaws                            |
| 14. | 8-011T     | Barraquer Needle Holder, Small                       | 33. | 12-325*         | Pick Vitreoretinal Forceps                                     |
|     | 8-031T*    | Barraquer Needle Holder, Medium                      | 34. | 13-084          | Membrane Spatula 2.5 mm  |
|     | 8-061T*    | Barraquer Needle Holder, Long Size                   |     | 13-092*         | Membrane Scratcher   |
| 15. | 8-025T     | Barraquer Needle Holder, Extra Fine Jaws, Small      | 35. | 12-5101         | Anterior Vitrectomy Cutter, B&L®                               |
|     | 8-045T*    | Barraquer Needle Holder, Extra Fine Jaws, Medium     |     | 12-5068*        | Anterior Vitrectomy Cutter, Alcon®                             |
|     | 8-051T*    | Barraquer Needle Holder, Extra Fine Jaws, Long       | 36. | 12-6000         | Titanium Backflush Actuator                                    |
| 16. | 11-040S    | Westcott Tenotomy Scissors, Blunt Tips, 13 mm Blades | 37. | 12-5017         | Brush Tip Cannula  |
|     | 11-042S*   | Westcott Tenotomy Scissors, Blunt Tips, 6 mm Blades  |     | 12-5063*        | Soft Tip Cannula   |
| 17. | 11-044S    | Westcott Stitch Scissors, Sharp Tips, 6 mm Blades    |     | 12-5151*        | Charles Flute Cannula  |
|     | 11-046S*   | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    |     | 12-5193*        | Diamond Dusted Soft Tip Cannula                                |
|     | 11-047S*   | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    | 38. | 12-5116         | Self-Retaining Infusion Cannula                                |
| 18. | 11-100S    | Knapp Straight Strabismus Scissors                   | 39. | 12-5191         | Dual Bore PFC Cannula  |
|     | 11-101S*   | Knapp Curved Strabismus Scissors                     | 40. | 12-5072         | Fiber Optic Probe, Alcon®                                      |
| 19. | 14-040T    | Lieberman Temporal Speculum                          |     | 12-5076*        | Fiber Optic Probe, B&L®  |
| 20. | 14-022S    | Barraquer Wire Speculum, Adult Size                  | 41. | SmartSil5000*** | Purified Silicone Oil 5000 cSt                                 |
|     | 14-023S*   | Barraquer Wire Speculum, Child Size                  |     |                 |  |
| 21. | 15-301/303 | Silicone Bulb With Adapter                           |     |                 |  |
| 22. | 16-081S    | Towel Forceps  |     |                 |  |


ALBA PART PRINTER  
0300-4151546



- |     |            |  |     |                 |  |
|-----|------------|--|-----|-----------------|--|
| 1.  | 1-020S     | Bard Parker Handle                                   | 22. | 16-081S         | Towel Forceps  |
| 2.  | 1-010T     | Blade Holder   | 23. | 16-111S         | Schocket Double Ended Scleral Depressor                        |
| 3.  | 2-010T     | Castroviejo Caliper                                  | 24. | 18-304*         | Plastic Sterilization Tray, Single Level, Large                |
| 4.  | 2-100S     | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 25. | 18-305*         | Plastic Sterilization Tray, Double Level, Large                |
|     | 2-100T*    | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 26. | 18-307*         | Mini Sterilization Tray with Silicone Mat, 38 x 68,6 x 25,4 mm |
|     | 2-101T*    | Braunstein Fixed Caliper, 3.0-3.5 mm                 | 27. | 12-5173-1       | Reusable Two Step Trocar System with closure valves            |
| 5.  | 4-0600T    | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth       | 28. | 12-001T**       | Vitreoretinal Instrument Handle, One Finger Control            |
| 6.  | 4-0601T    | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth        |     | 12-003T**       | Vitreoretinal Instrument Handle, Squeeze Model                 |
| 7.  | 4-0607S    | Bishop-Harmon Suturing Forceps, 3 mm 1x2 Teeth       | 29. | 12-202-23       | Vertical Vitreoretinal Scissors, 70 Degree                     |
| 8.  | 4-091S     | Kelman-McPherson Tying Forceps, 8 mm Platforms       | 30. | 12-209-23       | Curved Subretinal Scissors, Curvature Radius 12 mm             |
|     | 4-092T*    | Kelman-McPherson Tying Forceps, 10 mm Platforms      | 31. | 12-301-23       | Diamondized End-Gripping Vitreoretinal Forceps                 |
|     | 4-174T*    | McPherson Angled Tying Forceps, 8 mm Platforms       |     | 12-4012*        | Vitreoretinal End Gripping Forceps With Micro Jaws             |
|     | 4-175T*    | McPherson Angled Tying Forceps, 10 mm Platforms      |     | 12-402-23*      | Vitreoretinal ILM «Nail» End Gripping Forceps                  |
| 9.  | 4-2201T    | Watzke Sleeve Spreading Forceps                      | 32. | 12-410-23       | Eckardt End-Gripping Forceps                                   |
| 10. | 4-2300T    | Bonaccolto Utility Forceps                           |     | 12-411-23*      | Tano Asymmetrical End-Gripping Forceps                         |
| 11. | 4-122S     | Halsted Mosquito Forceps, Straight, Long             |     | 12-420-23*      | End Grasping Forceps   |
|     | 4-123S*    | Halsted Mosquito Forceps, Curved, Long               | 33. | 12-321-23       | Spring Gripping Forceps  |
| 12. | 5-042      | Graefe Muscle Hook, Size 2                           | 34. | 12-325-23       | Pick Vitreoretinal Forceps                                     |
| 13. | 5-060      | Gass Retinal Detachment Hook                         | 35. | 13-097-23       | Delicate Membrane Pick   |
| 14. | 8-011T     | Barraquer Needle Holder, Small                       | 36. | 12-5181         | Disposable Two Step Trocar System                              |
|     | 8-031T*    | Barraquer Needle Holder, Medium                      | 37. | 12-5171         | Vitrectomy Cutter, B&L®  |
|     | 8-061T*    | Barraquer Needle Holder, Long Size                   |     | 12-5172*        | Vitrectomy Cutter HIGH SPEED, Alcon®                           |
| 15. | 8-025T     | Barraquer Needle Holder, Extra Fine Jaws, Small      | 38. | 12-6000         | Titanium Backflush Active Aspiration                           |
|     | 8-045T*    | Barraquer Needle Holder, Extra Fine Jaws, Medium     | 39. | 12-5162         | Brush Tip Cannula  |
|     | 8-051T*    | Barraquer Needle Holder, Extra Fine Jaws, Long       |     | 12-5161*        | Soft Tip Cannula   |
| 16. | 11-040S    | Westcott Tenotomy Scissors, Blunt Tips, 13 mm Blades |     | 12-5164*        | Charles Flute Cannula  |
|     | 11-042S*   | Westcott Tenotomy Scissors, Blunt Tips, 6 mm Blades  |     | 12-5192*        | Diamond Dusted Soft Tip Cannula                                |
| 17. | 11-044S    | Westcott Stitch Scissors, Sharp Tips, 6 mm Blades    | 40. | 12-5144         | Self-Retaining Infusion Cannula, 500 µm                        |
|     | 11-046S*   | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    |     | 12-5203         | Dual Bore PFC Cannula  |
|     | 11-047S*   | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    | 41. | 12-5203         | Fiber Optic Probe, Alcon®                                      |
| 18. | 11-100S    | Knapp Straight Strabismus Scissors                   | 42. | 12-5177         | Fiber Optic Probe, B&L®  |
|     | 11-101S*   | Knapp Curved Strabismus Scissors                     |     | 12-5179*        |  |
| 19. | 14-040T    | Lieberman Temporal Speculum                          | 43. | SmartSil5000*** | Purified Silicone Oil 5000 cSt                                 |
| 20. | 14-022S    | Barraquer Wire Speculum, Adult Size                  |     |                 |  |
|     | 14-023S*   | Barraquer Wire Speculum, Child Size                  |     |                 |  |
| 21. | 15-301/303 | Silicone Bulb With Adapter                           |     |                 |  |

ALBARA PRINTING  
0300-4151546





- |     |          |  |     |                 |  |
|-----|----------|--|-----|-----------------|--|
| 1.  | 1-020S   | Bard Parker Handle                                   | 19. | 14-040T         | Lieberman Temporal Speculum  |
| 2.  | 1-010T   | Blade Holder   | 20. | 14-022S         | Barraquer Wire Speculum, Adult Size                                      |
| 3.  | 2-010T   | Castroviejo Caliper                                  |     | 14-023S*        | Barraquer Wire Speculum, Child Size                                      |
| 4.  | 2-100S   | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 21. | 15-301/303      | Silicone Bulb With Adapter   |
|     | 2-100T*  | Braunstein Fixed Caliper, 3.5-4.0 mm                 | 22. | 16-081S         | Towel Forceps  |
|     | 2-101T*  | Braunstein Fixed Caliper, 3.0-3.5 mm                 | 23. | 16-111S         | Schocket Double Ended Scleral Depressor                                  |
| 5.  | 4-0600T  | Castroviejo Suturing Forceps, .12 mm 1x2 Teeth       | 24. | 18-304*         | Plastic Sterilizing Tray, Single Level, Large                            |
| 6.  | 4-0601T  | Castroviejo Suturing Forceps, .3 mm 1x2 Teeth        | 25. | 18-305*         | Plastic Sterilizing Tray, Double Level, Large                            |
| 7.  | 4-0607S  | Bishop-Harmon Suturing Forceps, 3 mm 1x2 Teeth       | 26. | 18-307*         | Mini Sterilization Tray with Silicone Mat, 38 x 68,6 x 25,4 mm           |
| 8.  | 4-091S   | Kelman-McPherson Tying Forceps, 8 mm Platforms       | 27. | 12-001T**       | Vitreoretinal Instrument Handle, One Finger Control                      |
|     | 4-092T*  | Kelman-McPherson Tying Forceps, 10 mm Platforms      |     | 12-003T**       | Vitreoretinal Instrument Handle, Squeeze Model                           |
|     | 4-174T*  | McPherson Angled Tying Forceps, 8 mm Platforms       | 28. | 12-2029         | Vertical Vitreoretinal Scissors, 45 Degree                               |
|     | 4-175T*  | McPherson Angled Tying Forceps, 10 mm Platforms      | 29. | 12-2099         | Curved Vitreoretinal Scissors  |
| 9.  | 4-2201T  | Watzke Sleeve Spreading Forceps                      | 30. | 12-3019         | Diamondized Straight Gripping Vitreoretinal Forceps                      |
| 10. | 4-2300T  | Bonaccolto Utility Forceps                           | 31. | 12-4089         | Fine End-Gripping Forceps  |
| 11. | 4-122S   | Halsted Mosquito Forceps, Straight, Long             |     | 12-411-25*      | Tano Asymmetrical End-Gripping Forceps                                   |
|     | 4-123S*  | Halsted Mosquito Forceps, Curved, Long               |     | 12-420-25*      | End Grasping Forceps   |
| 12. | 5-042    | Graefe Muscle Hook, Size 2                           | 32. | 12-3259         | Pick Vitreoretinal Forceps   |
| 13. | 5-060    | Gass Retinal Detachment Hook                         | 33. | 13-0979         | Delicate Membrane Pick   |
| 14. | 8-011T   | Barraquer Needle Holder, Small                       | 34. | 12-5150         | Trocar System  |
|     | 8-031T*  | Barraquer Needle Holder, Medium                      | 35. | 12-5140         | Vitrectomy Cutter, 25psi - B&L®, DORC® Harmony, Nidek® CV-12000, VT-5000 |
|     | 8-061T*  | Barraquer Needle Holder, Long Size                   |     | 12-5157*        | Accurus HIGH SPEED Vitrectomy Cutter                                     |
| 15. | 8-025T   | Barraquer Needle Holder, Extra Fine Jaws, Small      | 36. | 12-6000         | Alcon®: Accurus  |
|     | 8-045T*  | Barraquer Needle Holder, Extra Fine Jaws, Medium     | 37. | 12-5160         | Titanium Backflush Active Aspiration                                     |
|     | 8-051T*  | Barraquer Needle Holder, Extra Fine Jaws, Long       |     | 12-5152*        | Brush Tip Cannula  |
| 16. | 11-040S  | Westcott Tenotomy Scissors, Blunt Tips, 13 mm Blades |     | 12-5156*        | Soft Tip Cannula   |
|     | 11-042S* | Westcott Tenotomy Scissors, Blunt Tips, 6 mm Blades  | 38. | SmartSil5000*** | Charles Flute Cannula  |
| 17. | 11-044S  | Westcott Stitch Scissors, Sharp Tips, 6 mm Blades    |     |                 | Purified Silicone Oil 5000 cSt   |
|     | 11-046S* | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    |     |                 |  |
|     | 11-047S* | Westcott Stitch Scissors, Sharp Tips, 3 mm Blades    |     |                 |  |
| 18. | 11-100S  | Knapp Straight Strabismus Scissors                   |     |                 |  |
|     | 11-101S* | Knapp Curved Strabismus Scissors                     |     |                 |  |

ALBARIARTPRINTER  
0300-4151546